


EIC DIRC

Greg Kalicy




March 23, 2015

Outline

- **Basics**
 - DIRC concept
 - Reconstruction approach
- **Main Goals**
 - Investigate possibility of pushing state-of-the-art performance
 - Feasibility of using DIRC in EIC detector
 - Integration of DIRC with other systems
- **Ongoing activities**
 - Simulations with DrcProp and Geant
 - Prototype in particle beam
 - Component lab tests


Synthetic fused silica prototype bar



Outline

- Basics
 - DIRC concept
 - Reconstruction approach
- Main Goals
 - Investigate possibility of pushing state-of-the-art performance
 - Feasibility of using DIRC in EIC detector
 - Integration of DIRC with other systems


- Ongoing activities
 - Simulations with DrcProp and Geant
 - Prototype in particle beam
 - Component lab tests



Outline

- **Basics**
 - DIRC concept
 - Reconstruction approach
- **Main Goals**
 - Investigate possibility of pushing state-of-the-art performance
 - Feasibility of using DIRC in EIC detector
 - Integration of DIRC with other systems
- **Ongoing activities**
 - Simulations with DrcProp and Geant
 - Prototype in particle beam
 - Component lab tests

Simulation



Outline


- **Basics**
 - DIRC concept
 - Reconstruction approach
- **Main Goals**
 - Investigate possibility of pushing state-of-the-art performance
 - Feasibility of using DIRC in EIC detector
 - Integration of DIRC with other systems
- **Ongoing activities**
 - Simulations with DrcProp and Geant
 - **Prototype in particle beam**
 - Component lab tests

Full system prototype in CERN test beam campaign



Outline


- **Basics**
 - DIRC concept
 - Reconstruction approach
- **Main Goals**
 - Investigate possibility of pushing state-of-the-art performance
 - Feasibility of using DIRC in EIC detector
 - Integration of DIRC with other systems
- **Ongoing activities**
 - Simulations with DrcProp and Geant
 - Prototype in particle beam
 - Component lab tests



Outline



- **Basics**
 - DIRC concept
 - Reconstruction approach
- **Main Goals**
 - Investigate possibility of pushing state-of-the-art performance
 - Feasibility of using DIRC in EIC detector
 - Integration of DIRC with other systems
- **Ongoing activities**
 - Simulations with DrcProp and Geant
 - Prototype in particle beam
 - Component lab tests

Setup to test photosensors in high magnetic field



Cherenkov Detectors


- Main Cherenkov detector concepts in particle physics:




- Compare ring image with expected image for $e/\mu/\pi/K/p$ (likelihood test) or calculate mass from track β using independent momentum measurement (B field, tracking).

DIRC Detection of Internally Reflected Cherenkov Light


- Charged particle traversing radiator (with refractive index n) with velocity $\beta = \frac{v}{c}$ emits Cherenkov photons on cone with half opening angle: $\cos\theta = \frac{1}{\beta n}$
- Cherenkov angle conserved during internal reflections of propagating photons.
- Photons exit radiator bars through focusing elements into expansion volume and are imaged on photon detector array.
- Photon detector array measures **x, y and time** of photons that exit radiator and defines **θ_c , φ_c and time of propagation of individual Cherenkov photons.**



Synthetic fused silica prototype bar





DIRC Example hit pattern





→ DIRC hit patterns do not look like **typical RICH detector**.

DIRC Example hit pattern




- DIRC hit patterns do not look like **typical RICH detector**.
- Part of the ring escapes, not totally internally reflected.
- Ring image gets folded due to propagation in bar/plate.


DIRC Example hit pattern



- DIRC hit patterns do not look like typical RICH detector.
- Part of the ring escapes, not totally internally reflected.
- Ring image gets folded due to propagation in bar/plate.
- Additional folding in the prism expansion volume.




DIRC Occupancy plot from experiment





122° polar angle

- Comparison of the **pixelized simulated data** to **test beam data**.
- Grey dots in the background are true hit positions from the simulation.





DIRC Example hit pattern

- What is so challenging here?
 - Ring segments corresponding to different particles are close and overlapping




DIRC Example hit pattern

- What is so challenging here?
 - Ring segments corresponding to different particles are close and overlapping





Pions vs **Protons** at 3 GeV/c.
210k accumulated events.
(Hit pattern from simulation)




DIRC Example hit pattern

- What is so challenging here?
 - Ring segments corresponding to different particles are close and overlapping
 - Pixelated image (+ additional background)




Pions vs Protons at 3 GeV/c.
210k accumulated events
(Pixelated Hit pattern from simulation)




DIRC Example hit pattern

- What is so challenging here?
 - Ring segments corresponding to different particles are close and overlapping
 - Pixelated image (+ additional background)
 - Reconstruction from 20-50 photons per event



Pions at 10 GeV/c. Example of 6 single events.
(Pixelated Hit pattern from simulation)



Barrel DIRC Expected performance




- Difference in Cherenkov angle between π and K at 4 GeV/c is 6.5 mrad
- Required per track resolution has to be 2 mrad or better

$$\sigma_{\Theta_C}^{track} = \sqrt{\left(\sigma^{corelated}\right)^2 + \left(\frac{\sigma_{\Theta_C}^{photon}}{\sqrt{N_{pe}}}\right)^2}$$

- $\sigma^{corelated} < 1.5 - 2 \text{ mrad}$ **Correlated term:**
tracking detectors, multiple scattering
- $\sigma_{\Theta_C}^{photon} < 8 - 9 \text{ mrad}$ **Single photon Cherenkov angle resolution:**
bar size, pixel size, chromatic, bar imperfections
- $N_{pe} > 20$ **Number of photons:**
bar size, bar imperfections, photon detection efficiency of the detector




DIRC Cherenkov angle reconstruction method

- Reconstruction method**
 Pixel position + bar location define photon direction at bar end, stored in Look-up table, combined with particle track to calculate Θ_C .
- Photon path not unique**
 Additional combinatorial background in Θ_C



DIRC Cherenkov angle reconstruction method

- **Reconstruction method**
Pixel position + bar location define photon direction at bar end, stored in Look-up table, combined with particle track to calculate Θ_C .
- **Photon path not unique**
Additional combinatorial background in Θ_C



Improving DIRC Performance

- **Make DIRC less sensitive to background**
 - decrease size of expansion volume
 - use photon detectors with smaller pixels and faster timing
 - place photon detector inside magnetic field
- **Investigate alternative radiator shapes** (plates, disks)
- **Push DIRC π/K separation by improving single-photon θ_c resolution**

BABAR-DIRC Cherenkov angle res.: **9.6 mrad** per photon → **2.4 mrad/track**

Limited in BABAR by:

- size of bar image
- size of PMT pixel
- chromaticity ($n=n(\lambda)$)

~4.1 mrad
~5.5 mrad
~5.4 mrad


Could be improved for new DIRCs via:

- focusing optics
- smaller pixel size
- better time resolution


9.6 mrad → **4-5 mrad** per photon → **< 1.5–2 mrad/tr.**

DIRC Design Options


Expansion volume: oil tank/prism



Radiator: narrow bar/wide plate




Focusing system: different lenses




Monte Carlo Studies

- **DrcProp: stand-alone package for ray tracing simulations includes:**
 - Detector geometries
 - Beam properties
- **Stand alone Geant4 simulation package**
 - Physical processes
 - Will be integrated with MEIC
- **MEIC simulation**
 - EIC environment

Example of event display from DrcProp






Geant4 simulation of plate and prism geometry






Focusing Prototype of 3-component lens

- Prototype of 3 component high refractive lens without air gap.
- Lens produced and already tested in particle beam with PANDA Barrel DIRC group




Focusing Simulation of 3-component lens



Procured 3-component lens, Planacon MCP
Procured 3-component lens, 2x2 mm MCP




- A high-performance DIRC will require new small-pixel photosensors.
- The simulation confirms that a resolution close to 1 mrad (6 GeV/c) can be reached at forward (and backward) angles.

To be continued...



Other DIRCs BELLE II TOP


- **Belle II TOP** is barrel DIRC-type RICH with emphasis on fast timing
- **PID goal:** 3σ π/K separation for $p < 4$ GeV/c
- **Radiator:** fused silica plate 45cm wide, 2cm thick, 250cm long.
- TOP barrel formed by 16 plates
- Small expansion volume (10cm depth)
- **Photon detector:** array of 32 Hamamatsu SL-10 MCP-PMTs per sector, 512 in total




Barrel DIRCs



**BABAR
DIRC**





**BELLE II
TOP**



**PANDA
BARREL
DIRC**

Radiator geometry	Narrow bars (35mm)	Wide plates (450mm)	Narrow bars (32mm)
Barrel radius	85cm	115cm	48cm
Bar length	490cm (4×122.5cm)	250cm (2×125cm)	240cm (2×120cm)
Number of long bars	144 (12×12 bars)	16 (16×1 plates)	80 (16×5 bars)
Expansion volume	110cm, ultrapure water	10cm, fused silica	30cm, mineral oil
Focusing	None (pinhole)	Mirror	Lens system
Photon detector	~11k PMTs	~8k MCP-PMT pixels	~15k MCP-PMT pixels
Timing resolution	~1.7ns	<0.1ns	~0.1ns
Pixel size	25mm diameter	5.5mm×5.5mm	6.5mm×6.5mm
PID goal	3 s.d. π/K to 4 GeV/c	3 s.d. π/K to 4 GeV/c	3 s.d. π/K to 3.5 GeV/c
Timeline	1999 - 2008	Installation 2015	Installation 2017/18

DIRC Example of number of hits per track




~20 photons per track predicted by simulation at 124° polar angle, 10 GeV/c momentum beam.

- 908 photons generated per track
- ~59% of photons propagate to the end of the bar.
- ~25% enter expansion volume.
- ~19% reach photo cathode of MCP-PMTs.
- ~2.2% detected.

DIRC Focusing options

128° Polar Angle

90° Polar Angle




- No focusing, bar directly coupled to the prism

- Standard UV coated lens with 2 mm air gap between bar and prism.


- Compound lens Fused silica/NLaK without air gap.

a) No Focusing




Test beam data

b) Standard lens



Test beam data

c) Compound lens




Test beam data

DIRC Focusing options

128° Polar Angle

90° Polar Angle



- In this configuration 10% improvement of photon yield, consistent with simulations.
- For particle polar angle close to 90° difference much more dramatic.

Number of hits per track:

