Thermal Neutron Cross Section Measurements

NDAG Meeting November 19, 2021

D. Fritz, Y. Danon

This work was supported by the Nuclear Criticality Safety Program, funded and managed by the National Nuclear Security Administration for the Department of Energy.







Overview

• <u>Motivation</u>:

- 1. Many materials lack high quality total cross section measurements from 0.5 meV 1 eV required for validation of thermal scattering library (TSL) evaluations.
- 2. Evaluations in this region need immediate validation as thermal scattering can dramatically alter the neutron multiplication factor of a system.



Completed Work:

- 1. Designed and developed cold moderator capability for use in moderating neutrons to the subthermal region (<20 meV).
- 2. Performed thermal cross section measurements that encompass the entire thermal region (0.5 meV 1 eV) for polyethylene, polystyrene, Plexiglas G/G-UVT and yttrium hydride.







Design Specifications

- Neutronics:
 - Polyethylene
 - Large enough to cover collimation
 - As close to source of neutrons as possible
 - As cold as possible in a consistent fashion
 - Maximize gain over ETT alone (existing target)
 - Metric: Integral counts in 1-5 meV range
- Thermal/System:
 - HDPE
 - < 30 K w/ 1.5 kW on target
 - Good connection to cryostat
 - System
 - Must couple to ETT target
 - Long cold finger \rightarrow protect cold head
 - Large cold finger cross-sectional area
 - Portable; Minimal material; Can withstand Hi-Vac
 - Single Stage Cryo-cooler











Cryostat rolls in front of ETT







Design Choice



Design Gain Comparison (Over ETT Alone)



- Both Design I and Design II provide significant gain over the ETT.
- Despite the improvement of Design II over Design I, Design I was proceeded with due to safety considerations.







System Experimental Gain Validation



- After LINAC power normalization, the ETTC achieves a gain up to 8 over the ETT.
- ETTC extends flux below 1 meV.







Thermal Cross Section Validation



- Excellent agreement between ETT and ETTC cross section for polyethylene (PE).
- Good agreement between all RPI measured cross sections, other measured cross sections, and TSL evaluation.
- Experimental error bars account for all experimental sources of uncertainty and were calculated via a covariance matrix.







Polyethylene – Closer Look



*Excellent agreement between all RPI experiments







Slide 9 of 16

Polyethylene – Averaged Cross Section



Weighted average of all polyethylene measurements performed



•



• ENDF/B-VIII.0 evaluation is higher than measured cross section at most energies



Thermal Cross Section Measurements: Polymers



- Generally good agreement between measured cross section and TSL evaluation for polymers.
- Differences in cross section between Plexiglas G and Plexiglas G-UVT.
- Experimental error bars account for all experimental sources of uncertainty and were calculated via a covariance matrix.







Thermal Cross Section Measurements: YHx



- Generally good agreement between yttrium hydride measured cross section and TSL evaluation for both concentrations.
- Experimental error bars account for all experimental sources of uncertainty and were calculated via a covariance matrix.
- Estimated uncertainty on H/Y ratio from samples not shown here.







YHx – Closer Look



- Misalignment of Bragg edges between experiment and yttrium evaluation
- Cross-over in cross section seen in YH_{1.68} samples





- Potential explanations for differences:
 - Numerical difficulties in NJOY



Be Thermal Cross Section Measurement



- Discrepancies between measured Be cross section and evaluated data
- Discrepancies between NNDC ENDF/B-VIII.0 Phonon Spectrum and ENDF/B-VIII.0 ACE file
- Both NNDC and ENDF ACE had ENDF/B-VIII.0 Be capture data added.







Be – Closer Look



 Large differences in amplitude between evaluations and measurement is likely due to the crystalline structure of the Be samples used → evaluations assume no extinction for grain size in Be







Milestones

• FY19-21 Completed:

- Successful neutronic/thermal modeling and simulation of cold moderator system.
- Cryogenic system designed and procured.
- Tested cryogenic system steady state temperature of 29K with LINAC Power.
- Completion of cryogenic system characterization.
- Completion of total thermal cross section measurements for polyethylene, polystyrene, Plexiglas, yttrium hydride and beryllium.
- FY22:
 - 12/2021: Finalization of all data reduction for all measurements.
 - 2/2022: Completion of necessary journal publications.





