TPOT Gas System


Robert Pisani


January 26, 2022

Gas System Technical Overview

TPOT Gas System Requirements

- •Argon + 5% isobutane : mixture for good gain and stability, slow gas against Lorentz angle
- •Module Volume 1.4 liters
- •Total Volume 11.2 liters
- •Few turns per hour
 - •50ccm per module (3L/hr)
 - •400ccm total
- •Slight overpressure above atmosphere (1mBar 0.40"W.C)
- •Compensate for changes in barometric pressure
- •High purity system
- •Single Pass
- •Remote monitoring
- •Alarmed





1/26/2022


The TPOT Gas System


•Has mixing capability for 3 gases (95.0% R134a / 4.5% Isobutane / 0.5% SF6)

- •Will Switch to Argon mixture (Ar-Iso) (95-5)
- •Calibrated mass flow controllers will mix gas and monitor ratio
- •Isobutane analyzer will monitor mixture.
- •Flowrate ~ 400ccm
- •0-2mBar (~0.5" w.c.) above atmosphere

to +/- 1/8 w.c."

- •Single Pass system
- •Supply pressure from Rack ~ 30" W.C.
- •Return Pressure –Atmosphere
- •Analyzers to monitor mixture and impurities
- •In line purifiers if needed
- •High purity system •<500 ppm O2 •<100 ppm H2O
- •PC Controlled
- •Fully alarmed






1/26/2022

Original TOF.W System Parameters


	TOF.W	For TPOT
Mixture	95%R134a + 4.5%Iso +0.5% SF6	95%R134a + 5%Isobutane
Compressor pressure	40-80 " H2O	
Supply pressure	30"+/-0.05 H20	30"+/-0.05 H20
Return pressure	1" +/-0.05 H2O	Vent to Atmosphere
TOF.W Diff Pressure	1.5 +/05"WC	0.4" +/05"WC
Recirculation flow	650-1000 ccm	
Flow through chambers	Up to 850 ccm	Need 400ccm
Purge flow	2.0 l/m	2.01pm
Make-up mixture flow	100-300 ccm	
Oxygen content	< 500ppm	< 500ppm
Water content	<100ppm	<100ppm

TPOT Gas System Features

- Previously Operating PHENIX Gas system
- Basic Features
 - 95% Argon + 5% Isobutane
 - Mass flow controllers (calibrated every year)
 - Constant controlled pressure and flow
 - Straight through or recirculation
- Gas Monitoring
 - Monitors system and chamber pressures,
 - Flammable gas content, Independent Iso monitor
 - Monitor flow to chambers
 - O₂ and H₂O levels monitored
- PC controlled and monitored system
 - Collects and logs system parameters
 - Warnings and Alarms
 - Example: if Argon runs low, Iso shut off
 - Remote access



Major Components



CAI GPA2000 Gas purity analyzer with pressure compensation. Calibrated for Argon-Iso 0.1% accuracy

This will continually monitor the concentration of isobutane in our mixture.





Cermet II Hygrometer Measurement Range -0.001 to 9999 PPM



Teledyne InstaTrans XD O2 Monitor Range: 0-1000 PPM [configurable]




TESCOM™ ER5000 Series Electropneumatic Actuator---The ER5000 Series is a microprocessor based PID (Proportional, Integral, Derivative). Maintains return line pressure to better than 1/16" W.C.




The 300 Vue Mass FlowController. Precise control for accurate mixing. (0.2% of Full Scale) A high flow and low flow controllers are being used to get better accuracy in recirculation.

TPOT gas lines to IR




Existing Items




•Wiring exists

Some minor work to modify for TPot
All wiring and electrical components in Gas house are still intact.

•Some new runs to IR will be needed.



1/26/2022



Gas Summary

- PHENIX hosted many high purity gas systems (HBD detector, RPC, TOF, DC, PC, RICH, Mutr, Muid).
- The TPOT gas system will be the PHENIX TOF.W gas system which successfully ran for many years in PHENIX
- All major components exist.
- The system will have gas and contamination monitoring..
- For Flammables,
 - All mass flow controllers will be calibrated every year.
 - Calibrated flow controllers will be used to mix gas
 - A separate Isobutane analyzer will be run to monitor mixture
- The system will be fully alarmed with the ability of remote monitoring.