Field response check for APA1 w


Xuyang Ning & Wenqiang Gu 02/20/2025




Recap

- ☐ For APA1 w plane:
- ☐ Template "fit" the shape of the data by renormalizing the electron tracks.
 - ➤ "create" the shape by renormalizing the electron track at 0.4 and 0.5 pitch.
 - ➤ Add time stretch.


☐ Signal processing may not work well.







Simulation w/ and w/o CNR


Coherent noise removal not applied






Coherent noise removal applied







Data w/ and w/o CNR





Coherent noise removal not applied















Redo everything with no CNR




- 5 Parameters
 - Different renormalizing factors to 0.4 pitch (f_5) and 0.5 pitch (f_6)
 - One starting point for the 2 electron paths
 - Two time stretch factors for 2 electron paths

- A python script to control everything:
 - Minimize T_{total}^2
 - MCMC

Generate renormalized field response from wire-cell-python

Simulation using a wire-cell-tookit



Compare data with simulation in a root script:


Compute a test statistic T_{total}^2


Do fitting using MCMC

- Define T_{total}^2 :
 - For w plane, compare waveforms and power density spectra
 - For each plot:

$$T^{2} = \sum_{i=1}^{nbins} \frac{(x_{i} - m_{i})^{2}}{|m_{i}|} / nbins$$


- T_{total}^2 : For different plots in different tracks, multiply all the T^2
 - Only use 7 tracks, because the last one is strange





Do fitting using MCMC


- A very simple MCMC with the help of OpenAl:
 - Since the T_{total}^2 is dependent on simulation, the same parameters cannot produce same result all the time.
 - MCMC only relies on the result of the last step.
 - It is easy to deal with lots of parameters.
- Random walk: a new parameter set is randomly generated according to the current parameters.
 - The step size is fixed.
- ightharpoonup If $T_{total}^2(current\ parameters) > T_{total}^2(new\ parameters)$
 - Accept new parameters;
- > Else
 - > new parameters still has some chance to be accepted.
 - Metropolis-Hastings acceptance rule.





-0.04_2.10_1.00_2.50_650


Best one this time

waveform













Best one this time


waveform
















Raw waveform ANf APA1_w_28548_456790


raw_wf_ANF_1_w Entries 1.685634e+07






Best one this time

waveform







400

Raw waveform ANf APA1_w_28548_456790

raw_wf_ANF_1_w

Power spectra






200


300


100


-100


































Power spectra







150


200

Check decon result




Newly fitted field response







Check decon result



Newly fitted field response






Add a 'filter' to help ROI

- Overall_resp[1]
 - overall response that combine field response and electronic response.
- A 2D "filter" can be obtained by dividing these 2D response bin by bin.

https://github.com/WireCell/wire-cell-toolkit/blob/master/sigproc/src/OmnibusSigProc.cxx#L1054

```
for (size_t i = 0; i != overall_resp[plane].size(); i++) {
    for (int j = 0; j != m_fft_nticks; j++) {
        r_resp(i, j) = overall_resp[plane].at(i).at(j);
        // if(plane==1) {
        // std::cout<<r_resp(i, j)<<"\t";
        // }
    }
    // if(plane==1) {
        // std::cout<<std::endl;
        // }
}


// cxit(0);
// exit(0);</pre>
```



Add a 2D "filter"

- Best fit field response for charge decon
- Use w=0V field response for ROI finding
 - Apply the 2D "filter"
 - ➤ In fact we change the overall_resp[1] before ROI finding, and reload the original field response before calculate decon charge.


The ROI is wide but protect more region



Add a cut on decon charge

☐ Further add a threshold on decon charge

☐ I think this can be the final result for APA1 w plane field response.



- ☐ Tasks left:
 - > Find a way to include the "filter" into the configuration files.
 - Do some other relevant tests, like seeing the performance in wire-cell-imaging