





MPGD Barrel Outer Tracker (MPGD-BOT)

Kondo Gnanvo, Jefferson Lab – RD&I Group

ePIC MPGD Detector Subsystem Leader (MPGD DSL)

Incremental Design and Safety Review of the ePIC MPGD Trackers December 15th, 2025

JLab Team:

Xinzhan Bai, Seungjoon Lee on behalf of MPGD DSC

Charge Questions Addressed

- 1. Are the technical performance requirements appropriately defined and complete for this stage of the project?
- 2. Are the plans for the various sub-systems appropriately documented and complete for this stage of the project?
- 3. Are the current plans for the detector likely to achieve the technical performance requirements, with a low risk for cost increases, schedule delays, and technical problems?
- 4. Are the schedule assumptions for the fabrication of the various sub-systems and assembly plans reasonable and consistent with the overall detector schedule?
- 5. Have ESH&Q and QA considerations been adequately incorporated into the plans at their present stage?
- 6. Have the recommendations from previous reviews been adequately addressed

Outline

Roles & Requirements

Charge 1

MPGD-BOT Layout

Charge 1,3

Technology & Performances

Charge 3

❖ PED - Engineering Test Article


Charge 3, 4, 6

❖ Assembly Schedule, QA & ES&H

Charge 4, 5


Summary

Roles & Requirements




Primary role: Pattern recognition, additional hit point & redundancy

- * MPGD-BOT fast timing capability will help with the pattern recognition to discriminate between physics and background signals during track reconstruction.
 - Expected timing resolution \sim 10 ns → achievable with thin-gap hybrid MPGDs
- ❖ Will provide additional hit for robust track reconstruction and for redundancy
 - Expected spatial resolution:
 - ✓ \sim 70 µm for perpendicular tracks
 - ✓ ~150 µm for inclined tracks $[0, 45^{\circ}]$ → achievable with thin-gap hybrid MPGDs



Additional role: Enhance hpDIRC performance

- MPGD-BOT close to hpDIRC detector with acceptance matching
- Improved angular and space point resolution at the hpDIRC level






MPGD-BOT Layout


MPGD-BOT Layout

24 planar detector modules

- ❖ 2 sectors in z made of left and right planar modules
 - Keeping zones: z = [-164.5 cm, 174.5 cm]
- 12 modules in $(r \times \varphi)$ arranged in dodecagon geometry
 - Keeping zones: r = [73.5 cm, 76 cm]
 - No overlaps in r $\times \phi \rightarrow \sim 12 \%$ acceptance gap
- ❖ Total readout electronic channels ~86K





MPGD Barrel Outer Tracker (MPGD-BOT)





MPGD-BOT module


https://eic.jlab.org/Geometry/Detector/Detector-20240515102931.html











Service/Parts	Per Module	Total	Parameter	
Frontend Boards (FEBs)	14	336	256 Channel / FEB	
High Voltage	1 (through divider)	24	3.2 mm OD	
Low Voltage	14	336	6 mm OD	
Gas lines	2	48	3 mm ID, 4 mm OD	
Cooling	4	96	6 mm ID, 8 mm OD	
Data Cable	14	336	Optical Fiber	
Environment Sensor	2	48	Temp. & Humidity	
Ground			Depends on grounding plan	

- ❖ Parameters of service line are subject to change as FEB design progresses.
- ❖ The cooling requirement has not yet been finalized.
- ❖ The current cooling approach utilizes water cooling with a heatsink.
- ❖ MPGD-BOT, MPGD-ECT & CyMBaL will collaborate on cooling system design.

Technology & Performances


Thin-gap GEM-µRWELL hybrid detector

- ❖ Drift gap: 1.5 mm & transfer gap: 2.5 mm
- * Double amplification stage (GEM & μRWELL)
- ❖ Capacitive-sharing U-V strips readout layers
 - ❖ Pitch: 0.8 mm (3580 U & V strips / module)
- https://doi.org/10.1016/j.nima.2025.170791


Readout electronics:

- Front End Boards (FEBs) based on SALSA chips
- ❖ 14 FEBs / module (4 SALSA chips i.e 256 e-ch / FEB)


Cross section view of thin-gap GEM-µRWELL



CAD model of thin-gap GEM-µRWELL hybrid



Capacitive-sharing U-V strip readout – baseline design




https://doi.org/10.1016/j.nima.2022.167782

Test beam @ FNAL (2023) and @ JLab (2025) for performance study and optimization of thin-gap GEM-μRWELL prototypes


FNAL test beam (06/2023): Position resolution vs. track angle

- * 3 protos: 0.5 & 1 mm thin-gap GEM-μRWELL, 3-mm std μRWELL
- ❖ 4 HV scans for efficiency study with Ar:CO₂ (80:20) gas mixture
- Track angle scan $(0-45^{\circ})$ for position resolution comparison studies




JLab Hall D test beam (06/2025): Efficiency vs. gas mixtures

- * 2 protos: 1-mm and 1.5-mm thin-gap GEM-μRWELLs
- * HV scan for efficiency study various Ar-based mixtures
- ❖ Argon gas mixture: Ar:CO₂ & Ar:CO₂:iC₄H₁₀





EIC Generic Detector R&D Program: https://wiki.bnl.gov/eic/index.php?title=File:ERD_tgMPGD_FY22_endOfYearReport_final.pdf


- ❖ Spatial resolution improved by 2 w.r.t to standard gap MPGD.
- ❖ Efficiency: ~98% with proper choice of gas mixture
 - $Ar:CO_2:iC_4H_{10}$ (90:7:3) is the baseline choice.
- ❖ High degree of flexibility and versatility on HV settings
 - GEM, the μRWELL, drift and transfer field
 - Guaranteed operation stability at high performance




https://doi.org/10.1016/j.nima.2025.170791










PED: Engineering Test Article

MPGD-BOT module: Project Engineering Design Test Article


Charge 3




- ❖ All essential detector parts have been produced and are in hand at JLab
 - GEM & cathode foils
 - μRWELL + U-V strip readout PCB
 - Cathode and GEM Support frames
 - 30 Panasonic-to-Hirose connectors
- ❖ Next step: preparation for assembly
 - GEM HV sectors tests setup
 - GEM stretcher device
 - Assembly flow hood space
 - Cleaning & coating of the frames



μRWELL PCB front view



GEM foil



μRWELL PCB back view



Support frames

Preliminary Design and Safety Review of the EIC Tracking Detectors – Mar 20-21, 2024

MPGDs



Comment about charge 4

The production of the moderately high number of MPGD structures should not pose a major risk for the project. Still, the full-size structures (1800 mm long) should be produced as soon as possible, and the prototypes should be constructed to ensure timely realization of the prototyping and production plans.

Electron-Ion Collider

Incremental Preliminary Design and Safety Review of


Workflow for the Assembly

Components:


- Detector Frame
- GEM Foil
- μRWELL Foil

Includes cleanroom activity and non-cleanroom activity


Highly customized instruments – tailored to detector size and shape



MPGD test lab (EEL-126 @ JLab)




Clean room layout (EEL-122 @ JLab)



Both areas are tightly packed with equipment & instruments

GEM Foil HV Nitrogen Test Box


- HV test of unframed GEM
- HV test of Framed GEM
- HV test of the GEM after glued onto cathode





Stretch GEM Foil mechanical


- GEM glued to support frame
- Frame GEM glued to cathode / honeycomb





Detector assembly stand

- Glue framed GEM to µRWELL PCB
- Require a very clean space
- Laminar flow hood inside clean room





Electron-Ion Collider

PED Test Article: Frames - Sanding, Cleaning & Coating

Large ultrasonic tank for frame cleaning

- After sanding the frames
- Ultrasonic bath to clean the frame filled with de-ionized water
- Large volume ultrasonic thank
- Pending installation in EEL-126



Fume hood for frame coating

- Apply a thin film of varnish to protect the spacers after sanding and ultrasonic bath cleaning
- Varnish work requires a fume Hood
- Pending installation in EEL-126



Frame storage and dryer in clean room 121

- Dry the frames after USB cleaning
- Process needs 48 hours



Cutting the honeycomb to size





Inserting the honeycomb in the frame



Sandwich the honeycomb b/w G10 skins



Electron-Ion Collider

Assembly Plans, ES&H & QA

Module production sites:

- Jefferson Lab (JLab)
- University of Virginia (UVa)
- Florida Institute of Technology(FIT)

Vast experience with large scale and large MPGD projects among the 3 institutions

- ❖ 10 modules per site (8 + 2 spares)
- Overall effort under JLab supervision
- Common QA / QC procedures for all sites

Activities / schedule	01/2024 - 01/2025	01/2025 – 06/2026	03/2026 - 12/2026	01/2027 - 01/2028	01/2028 - 01/2031	01/31 -
PED Test article design	JLab					
PED Test article I		JLab				
PED Test article II			JLab			
Pre-production				JLab, UVa, FIT		
Production modules					JLab, UVa, FIT	
Commissioning @ BNL						JLab, UVa, FIT

All 3 assembly sites have fully equipped MPGD Detector Lab

- ❖ Fully equipped CLASS 1000 Clean rooms for module assembly
- * Cosmic tracking telescope setup with coincidence trigger counters and readout & DAQ system
- ❖ X-ray setup for high-rate studies and long-term stability ...
- ❖ Will setup DAQ and readout system for SALSA MPGD readout system


JLab MPGD Clean Room: New capacity for large MPGD module assembly



UVa Clean Room: SBS GEMs, MOLLER GEMs, PRad GEMs, CLAS12 μRWELL, Hall D GEM-TRD prototype



FIT Clean Room: assembly of CMS GE1/1 GE2/1 and ME0 GEMs



MPGD-BOT: Quality Assurance (QA) – ES&H

*** QA/QC**:

- HV test in N_2 of GEM and μ RWELL sectors before, during and after assembly (in clean room)
- Gas leak test and sealing after module assembly
- Electrical & capacitance test of the connectors: (QA at CERN), additional check at production sites
- Check of the pedestal noise and offset of all U & V strips after FEB connections
- Tag dead strips and replace FEB cards with dead readout channels
- Efficiency and relative gain uniformity studies with cosmic setup
- Absolute gain measurement with Fe55 (local measurement)

***** ES&H:

- Assembly will follow guidelines in compliance with ES&H requirements at JLab and the production sites
- Flammable gas Ar: CO_2 : iC_4H_{10} → Safety concern or implementation procedure will be implemented
- HV supply via simple passive divider → Procedure to isolate any electric point on the detector
- The ES&H procedures will be documented in the ePIC documentation database
- BNL ES&H procedures will be implemented during commissioning and installation in ePIC at BNL

Common gas mixture for all 3 MPGD subsystems: What does it take


Current gas mixture choice for the baseline detectors:

- **♦** MPGD-BOT baseline gas mixture: Ar:CO₂:iC₄H₁₀ (90:7:3)
 - Open to explore different ratio i.e.. 93:2:5 → will be studied in test beams at JLab & CERN in 2026
- ❖ CyMBaL baseline gas mixture: $Ar:CO_2:iC_4H_{10}$ (93:5:2) under consideration → likely not an option for MPGD-BOT
- ❖ MPGD-ECT is has been experimenting with Ar:CO₂:CF₄ (45:10:45)
 - CF_4 based gas mixture not under consideration for MPGD-BOT \rightarrow not required to achieve good timing


The Goal is to operate all three gaseous trackers with the same gas mixture

- We will agree on a mixture and ratio that satisfies all 3 trackers if achievable \rightarrow Ar:CO₂:iC₄H₁₀ mixture is a good starting point
- Several performance parameters need to be optimized simultaneously and sometimes they are contradictory
 - Fast gases (CF4-based) are **good for timing** but **bad for Lorentz effect** which also affects position & timing → see next slide
 - Large fraction of Ar is **good for gain** (critical for CyMBaL) but one must be careful about **stability** (i.e., increase discharge rate)
- ❖ More studies of different gas mixture with all 3 subdetectors will be performed in test beams in 2026 to converge to a common choice

Gas Selection: Garfield simulation of various gas mixtures



MPGD-BOT: Readout Electronics based on SALSA ASIC


















Short pigtail / on board

Fibers of adapted length between patch panels

Short pigtail / on board

FEB LV Powering scheme



- ❖ SALSA ASIC → same readout & DAQ as CyMBaL & MPGD-ECT
- ❖ Form factors of the on-detector FEB, DC-DC board, cooling scheme ... are subsystem specific
- ❖ JLab Electronics group will adapt FEB, DC-DC boards from CEA Saclay to MPGD-BOT specific needs


MPGD-BOT: Assembly plans & schedule

12/2025 01/2027 01/2028 01/2031 03/2031 Commissioning & PED & validation Production Pre-production Shipment to BNL Installation

- full size MPGD-BOT Test. article modules
- Procurement of GEM foils, μRWELL **06/2025**
- Assembly @ JLab &Test in beam @ CERN - 06/2026
- Have time for a second test article if needed before 2027
- Assembly of one preproduction module module#0 by each production site
- Setup of infrastructure and equipment in production sites
- Assembly and QA of 9 production modules at assembly sites
- Characterization with cosmic and radioactive sources
- Shipment of all 24 modules to BNL
- Commissioning at BNL with cosmic of all modules

- Commissioning
- Installation

High Level Installation Schedule From Elke's talk



Electron-Ion Collider

Incremental Preliminary Design and Safety Review of the ePIC N


31

Summary


- * ePIC MPGD-BOT Tracker will provide fast timing information for pattern recognition to support Si trackers in the barrel
- ❖ MPGD-BOT also provides hit redundancy capability and improved angular resolution to enhance hpDIRC performance
- * MPGD-BOT is based on thin-gap GEM-μRWELL hybrid detector technology
- ❖ Design of the full scale MPGD-BOT first test article module and procurement of essential parts are complete
- ❖ Setup of the clean room infrastructure for the assembly of the first test article at JLab is ongoing
- * Assembly is scheduled for the first quarter of 2026 and characterization in beam and with cosmic for summer 2026
- Plans and schedule for MPGD-BOT pre-production & production modules at the 3 assembly sites is scheduled for early 2027

Backup


MPGD-BOT module: Project Engineering Design Test Article




MPGD-BOT Parts: μRWELL foil design



MPGD-BOT Parts: GEM foil design



MPGD-BOT Parts: U-V strip readout PCB design




Technology choice: Thin-gap GEM-µRWELL Hybrid Detector

Challenges with standard (> 3-mm drift gap) MPGD

- Degradation of the spatial resolution with track angle.
- \bullet E × B in magnetic field negatively impact resolution

Development of Thin-gap MPGDs:


- ❖ Small drift gap improve spatial resolution at large angle
- \bullet Small gap \rightarrow minimize E × B effect in magnetic field
- Improve the detector timing performance



parametrization from EPJ Web of Conferences 174, 06005 (2018)

Thin-gap GEM-µRWELL detector concept

- hybrid amplification MPGD:
 - GEM (preamplification) and μ RWELL (main amplification)
 - Allow large detector gain and stable operating HV
- * Readout layer: 3-layer capacitive-sharing U-V strip readout
 - Achieve excellent spatial resolution with thin gap detector



cross-section view of thin-gap GEM-µRWELL detector


https://wiki.bnl.gov/eic/upload/ERD tgMPGD FY22 endOfYearReport final.pdf


Technology choice: Thin-gap GEM-µRWELL Hybrid Detector


Proof of concept

- Concept of thin-gap GEM-µRWELL
 hybrid prototype demonstrated in beam
 test at the Fermilab Test beam Facility in
 Summer 2023 (red plots)
- Space resolution < 150 μm and efficiency of 92% on average for 1-mm thin-gap GEM-μRWELL prototype (red dots) and for track in an angle range between 0 45 degrees.</p>
- Baseline technology for ePIC outerMPGD tracker


μRWELL + readout PCB









R&D funded by JLab administered DOE EIC Generic R&D Program as EICGENRandD_2022_23







https://wiki.bnl.gov/eic/upload/ERD tgMPGD FY22 endOfYearReport final.pdf

Design consideration: Some General considerations

Rate Capability

■ Not critical ~ 1 kHz/cm² or less

Radiation Hardness

- Not critical for the detectors
- Important for FEBs and RDO electronics boards
- ***** Temperature Stability
 - Not critical for the detector performances
 - Detector calibration should consider gas pressure variations
- Electronics power consumption and cooling
 - SALSA ASIC consumption ~ 15 mW/ channel at 1.2V
 - Air vs liquid cooling is under study at Saclay see Irakli's talk

