

DECal Calibration System

Analysis of the Spectrophotometric Flat Field System for the DECam Instrument in support of the Dark Energy Survey at the CTIO Blanco 4m Telescope

Presented by: William Wester¹ for the Dark Energy Survey Collaboration

in association with the DECal Team: Darren DePoy², Jennifer Marshall², Jean-Phillippe Rheault², Ting Li², Travis Prochaska², and David James³

¹Fermilab, ²Texas A&M University, ³Cerro Tololo Inter-American University


ABSTRACT

Data from the DECal spectrophotometric flat field system installed at the CTIO Blanco 4m telescope have been acquired throughout the first calendar year of operations of the DECam instrument. An initial analysis has been performed that allows for the measurement of the total system throughput (less atmosphere) – the response of the instrument as a function of wavelength, focal plane position, and time. The DECal hardware also includes an LED-based system for acquiring daily flat fields.

REQUIREMENTS


- Science requirements
 - Photometric calibrations: R10: spatially varying magnitude systematics < 2%, with a goal of < 1% for grizY and R11: < 0.5% on color zeropoints (1% for z-Y)
 R13: Synthesized grizY magnitudes agree with measurements to within 2%
- Specifications and Technical Requirements constrain optical corrector, filters, CCDs
 - TG12: Flat fielding procedure accurate to RMS < 0.5%
 - TO5: optical corrector optimized in performance approx. over 400-1000 nm
 TO8: exit pupil abosts < 3% gradient over 0.3 deg
 - TO8: exit pupil ghosts < 3% gradient over 0.3 deg
 TO15: Filter transmission > 85% in griz and TO16: < 10% peak-to-peak variation
 - TO20: Transmission shall fall within envelopes given in the Filter Specifications Document
 TD1: CCD nonlinearity < 1% with TD7: (g,r,i,z) CCD QE (60%, 75%, 75%, 65%)
- The DECal system is meant to help ensure and quantify the degree to which these specifications and requirements are met

HARDWARE




Hardware built by Texas A&M


LEDs for Daily Flat Fields





SYSTEM THROUGHPUT
















The total system throughput is the product of all the wavelength-dependent effects shown.

RAW DATA PRODUCTS


Typical scanning sequence

Point telescope at the flat field screen






Compute truncated averages in "super-pixels"




r-Filter response relative to none Jun/Feb double ratio demonstrating precision to better than 0.2%




Color represent one of eight different r-Filter scans obtained throughout the first DES year demonstrating repeatability and highlighting the rising and falling throughput edges and the out-of-bandpass light.

DERIVED DATA PRODUCTS

System throughput



Transformation equations



SOME ISSUES

Use of the DECal measurements requires care

- Illumination and Pupil Ghost effects are present. These are currently corrected in the DES data by using star flats
 A handful of amplifiers show some non-linearities at low light levels. In principle, each pixel could have it's own throughput ... what degree of segmentation is useful
- Full integration into a coadded image is complex with multiple epochs from overlapping exposures

CONCLUSIONS

Data from the DECal spectrophotometric flat field system are being collected and analyzed to verify technical specifications and requirements are being met at a level necessary for meeting the science requirements.

The system throughput is being measured and monitored as a function of wavelength, focal plane position, and time.