

# mRICH for EIC - past, present and future

Xiaochun He from Georgia State University on behalf of the EIC PID Consortium (eRD14)



# mRICH for EIC - past, present and future

Xiaochun He from Georgia State University

on behalf of the EIC PID Consortium (eRD14)

## Modular, Compact and Projective

#### **Abstract**



In this short presentation, I will briefly describe the mRICH design principle, its GEANT4 simulation and prototyping studies. A plan for future tasks will also be presented.

3/19/2020

#### **Abstract**



In this short presentation, I will briefly describe the mRICH design principle, its GEANT4 simulation and prototyping studies. A plan for future tasks will also be presented.

#### Following points will be addressed to the best of my knowledge

- Technology used: spell out clearly any risk associated, if any
- Momentum range covered: p versus theta and Nsigma vs. p
- Robustness of the design (e.g. sensitivity to magnetic field) and has a prototype been built?
- Are the electronics considerations clear (channel count, data size, rate, background)
- Time needed to complete the R&D and available workforce
- Status of Simulation and Reconstruction

#### **Abstract**




In this short presentation, I will briefly describe the mRICH design principle, its GEANT4 simulation and prototyping studies. A plan for future tasks will also be presented.

#### Following points will be addressed to the best of my knowledge

- Technology used: spell out clearly any risk associated, if any Aerogel + fresnel lens. NO tech risks!
- Momentum range covered: p versus theta and Nsigma vs. p 3 10 GeV/c k/pi, <2 GeV/c e/pi
- Robustness of the design (e.g. sensitivity to magnetic field) and has a prototype been built?
- Are the electronics considerations clear (channel count, data size, rate, background) YES
- Time needed to complete the R&D and available workforce YES, mainly to quantify the performance
- Status of Simulation and Reconstruction Well advanced GEANT4 simulation for the standalone mRICH








~ (aerogel thickness + lens focal length)

(Not to scale, for illustration purpose only)

## EIC mRICH – Working Principle





~ (aerogel thickness + lens focal length)

Aerogel Sensor lens

**Geant4 Simulation** 


With realistic material optical properties


(Not to scale, for illustration purpose only)

Prepared by X. He on 8/12/2019, EIC PID Consortium (eRD14 Collaboration)




Modular and compact ring imaging Cherenkov (mRICH) PID detector for EIC experiments








Modular and compact ring imaging Cherenkov (mRICH) PID detector for EIC experiments




New features: a) separation of optical and electronic components; b) longer focal length (6"); c) 3mm x 3mm photosensors.




2<sup>nd</sup> mRICH prototype was tested at Fermilab Test Beam Facility in June/July 2018









Modular and compact ring imaging Cherenkov (mRICH) PID detector for EIC experiments





- Projected e/pi separation of mRICH 2<sup>nd</sup> prototype detector (blue solid line)
- 2<sup>nd</sup> prototype detector can achieve 3sigma e/pi separation up to 2 GeV/c





- Projected K/pi separation of mRICH 2<sup>nd</sup> prototype detector (Green dots)
- 2<sup>nd</sup> prototype detector can achieve 3sigma K/pi separation up to 8 GeV/c

New features: a) separation of optical and electronic components; b) longer focal length (6"); c) 3mm x 3mm photosensors.




2<sup>nd</sup> mRICH prototype was tested at Fermilab Test Beam Facility in June/July 2018









Modular and compact ring imaging Cherenkov (mRICH) PID detector for EIC experiments







- Projected e/pi separation of mRICH 2<sup>nd</sup> prototype detector (blue solid line)
- 2<sup>nd</sup> prototype detector can achieve 3sigma e/pi separation up to 2 GeV/c




- Projected K/pi separation of mRICH 2<sup>nd</sup> prototype detector (Green dots)
- 2<sup>nd</sup> prototype detector can achieve 3-




New features: a) separation of optical and electronic components; b) longer focal length (6"); c) 3mm x 3mm photosensors.



2<sup>nd</sup> mRICH prototype was tested at Fermilab Test Beam Facility in June/July 2018







## **Beam Test Studies**


We had two beam tests already and are preparing for a third test at Lab in early May 2020.

### 1st mRICH Prototype Beam Test - Proof of Working Principle



Fermilab Beam Test Facility, April 2016






### 1st mRICH Prototype Beam Test - Proof of Working Principle



Fermilab Beam Test Facility, April 2016







C.P. Wong et. al. NIM A871 (2017) 13-19

#### 2nd mRICH Prototype - Improved Optical Component Design








- FLAT PANEL TYPE MULTIANODE PMT ASSEMBLY
  H13700 SERIES
- Longer focal length (Fresnel lens)
- 2. Smaller pixel size sensors

#### **FEATURES**


- High quantum efficiency: 33 % typ.
- High collection efficiency: 80 % typ.
- Single photon peaks detectable at every anode (pixel)
- Wide effective area: 48.5 mm × 48.5 mm
- 16 × 16 multianode, pixel size: 3 mm × 3 mm / anode




# The separation of the optical and electrical components in the improved mRICH design allows us to test different photosensors



Using four H13700 Multi-anode PMTs





Using three Hamamatsu SiPM Matrices







# The separation of the optical and electrical components in the improved mRICH design allows us to test different photosensors



#### Using four H13700 Multi-anode PMTs





To achieve the required PID separation power, the pixel size of photosensors should be 3mm x 3mm or smaller.



Each H13700 & SiPM matrix have 16 x 16 pixels (3mm x 3mm). Four sensors are needed to cover the imaging plane of mRICH. This leads to 1024 readout channels per module.

#### Using three Hamamatsu SiPM Matrices





## 2nd mRICH Beam Test - Verify the PID Capability



Fermilab Beam Test Facility, from July 25 to August 6, 2019






## 2nd mRICH Beam Test - Verify the PID Capability



Fermilab Beam Test Facility, from July 25 to August 6, 2019




120 GeV/c proton

## 2nd mRICH Beam Test - Verify the PID Capability




Fermilab Beam Test Facility, from July 25 to August 6, 2019





120 GeV/c proton



























## Position scans with 120 GeV/c proton beam





## Position scans with 120 GeV/c proton beam






## Outlook

## Construct more mRICH modules for performance test








- Preparing for a beam test at Jefferson Lab in early May of 2020 using secondary electrons.
- Plan to have another beam test at Fermilab once a tracking system is identified.
- Starting on mechanical design optimization.







- Modular
- Compact
- Projective

## mRICH modules in sPHENIX







**sPHENIX** 


#### **mRICH PID Parameterization**




• We are very close to implement the pid performance (see the plots on right) parameterization following Tom Hemmic's virtual PID base class: <mRICH.h>

```
#ifndef __MRICH_H__
#define __MRICH_H__
// Created on 3/15/2020 by Xiaochun He at Georgia State University
// for the EIC Yellow Report
// Wrapper class for mRICH (fast PID for the EIC mRICH)
#include "PID.h"
#include "mRICHPidFast.cxx"
class mRICH: public PID
public:
  mRICH(double trackResolution=0.5, double timePrecision=1.0, double pLow=3.0, double pHigh=10.0);
  virtual ~mRICH() {}
        valid (TVector3 p) {return (p.Mag() > pLow && p.Mag() < pHigh);}</pre>
  double numSigma(TVector3 p, PID::type PID);
                 (double numSigma, PID::type PID);
                 (double numSigma, PID::type PID) {return 3.0;}
                () {return myName;}
  void description ();
protected:
  std::string myName;
  mRICHPidFast pid:
  mRICHPidInfo info:
  double fTrackResolution; // resolution of the traker [mrad]
  double fTimePrecision: // time precision of the MCP-PMT [ns]
  double pLow;
  double pHigh;
  double fSensorQMefficiency; // photon sensor quantum efficiency
  int mRICH_ID: //
```

#### mRICH Performance based on GEANT4 simulation



- Projected e/pi separation of mRICH 2<sup>nd</sup> prototype detector (blue solid line)
- 2<sup>nd</sup> prototype detector can achieve 3sigma e/pi separation up to 2 GeV/c



- Projected K/pi separation of mRICH 2<sup>nd</sup> prototype detector (Green dots)
- 2<sup>nd</sup> prototype detector can achieve 3sigma K/pi separation up to 8 GeV/c

## **Summary**




- First of all, I would like to take this opportunity to thank the members of the EIC Detector R&D Committee for their continued support of the eRD14 PID Consortium on developing RICH-based PID detector technologies.
- Two successful beam tests of the mRICH prototypes demonstrated the mRICH working principles and validated the realistic implementation of mRICH in GEANT4 simulation.
- Further beam tests are needed in order to quantify the PID performance with improved detector design (including feasible photosensors and readout) and tracking capability.
- A array of mRICH modules have been implemented in the simulation framework of the sPHENIX experiment.
   THANKS!



# Backup slides


#### p/A h-endcap 10x100 GeV $Q^2 > 1 \text{ GeV}^2$ 10<sup>3</sup> barrel e-endcap $10^{2}$ **Dual-radiator** RICH 10 DIRC Aerogel RICH <sup>6</sup> rapidity 0 2




#### **PID Needs RICH Detectors**



**17** 







#### **PID Needs RICH Detectors**



- h-endcap: A RICH with two radiators (gas + aerogel) is needed for π/K separation up to ~50 GeV/c
- e-endcap: A compact modular aerogel RICH (mRICH)
   which can be projective π/K separation up to ~10 GeV/c
- barrel: A high-performance DIRC provides a compact and cost-effective way to cover the area.
   π/K separation up to ~6-7 GeV/c
- TOF (and/or dE/dx in TPC): can cover lower momenta.
- Photosensors and electronics: need to match the requirements of the new generation devices being developed both for the final systemand during the R&D phase







#### **PID Needs RICH Detectors**



- h-endcap: A RICH with two radiators (gas + aerogel) is needed for π/K separation up to ~50 GeV/c
- e-endcap: A compact modular aerogel RICH (mRICH)
   which can be projective π/K separation up to ~10 GeV/c
- barrel: A high-performance DIRC provides a compact and cost-effective way to cover the area.
   π/K separation up to ~6-7 GeV/c
- TOF (and/or dE/dx in TPC): can cover lower momenta.
- Photosensors and electronics: need to match the requirements of the new generation devices being developed – both for the final systemand during the R&D phase





#### **PID Needs RICH Detectors**




- h-endcap: A RICH with two radiators (gas + aerogel) is needed for π/K separation up to ~50 GeV/c
- e-endcap: A compact modular aerogel RICH (mRICH)
   which can be projective π/K separation up to ~10 GeV/c
- barrel: A high-performance DIRC provides a compact and cost-effective way to cover the area.
   π/K separation up to ~6-7 GeV/c
- TOF (and/or dE/dx in TPC): can cover lower momenta.
- Photosensors and electronics: need to match the requirements of the new generation devices being developed – both for the final systemand during the R&D phase

## mRICH – lens-based focusing aerogel detector design




#### Smaller, but thinner ring improves PID performance and reduces length




## mRICH – lens-based focusing aerogel detector design




#### Smaller, but thinner ring improves PID performance and reduces length



# Ring image from proton beam at an angle (110)






11 degree tilt downward

120 GeV/c proton beam



## mRICH readout with SiPM matrix sensors




#### Ring Radius and Number of Cherenkov Photons






No precision tracking was available. Beam size is ~6mm in radius.



(b) (c)







### **Ring Radius and Number of Cherenkov Photons**






No precision tracking was available. Beam size is ~6mm in radius.






(c)

#### Simulation







#### Ring Radius and Number of Cherenkov Photons






No precision tracking was available. Beam size is ~6mm in radius.







#### Simulation

