

mRICH for EIC YR

Xiaochun He, Murad Sarsour, Marco Contalbrigo and Zhiwen Zhao on behalf of the EIC PID Consortium (eRD14)



mRICH for EIC YR

Xiaochun He, Murad Sarsour, Marco Contalbrigo and Zhiwen Zhao on behalf of the EIC PID Consortium (eRD14)


YR progress report on May 1, 2020



Scope of the present mRICH parameterization for YR

To a first order estimation of the external requirement (i.e. tracking) for achieving the mRICH performance, our discussion will be mainly focused on a special case which will allow us to calculate the mRICH properties analytically.

This special case is that the particle is incident normal to the front of mRICH at the center. The mRICH performance with varying momentum direction and magnitude are part of the ongoing work which includes GEANT4 simulation and prototyping studies.






mRICH Internal Characteristics

Design tooling: analytical calculation, GEANT4 simulation and prototyping

See Murad Sarsour's presentation on fast parameterization



EIC mRICH – Working Principle




~ (aerogel thickness + lens focal length)

(Not to scale, for illustration purpose only)



EIC mRICH – Working Principle



~ (aerogel thickness + lens focal length)

Aerogel lens Sensor


Geant4 Simulation


With realistic material optical properties

(Not to scale, for illustration purpose only)

2nd mRICH Prototype (readout electronics not included)







- 1. Longer focal length (Fresnel lens)
- 2. Smaller pixel size sensors



FEATURES

- High quantum efficiency: 33 % typ.
- High collection efficiency: 80 % typ.
- Single photon peaks detectable at every anode (pixel)
- Wide effective area: 48.5 mm × 48.5 mm
- 16 × 16 multianode, pixel size: 3 mm × 3 mm / anode





mRICH Key Components - Optical (per module)

Components	Functions	Specs/Requirements	Risk	Mitigation
Mechanical (shoe box)	Supporting frame structure and light blocking.	Dimension: 11cm x 11cm x 20cm. Need mounting fixture designed for installing an array of mRICH modules.	None	Light and steady
Fresnel lens	Imaging focusing, UV filtering	Focal length ~6" (acrylic for UV filtering)	None	Good optical characterization scheme will be needed.
Mirrors	Reduce Cherenkov photon loss when particle is incident at a larger angle.	Rough dimension: Trapezoid shape: 15cm (height) x 11.5cm (width1) x 10cm (width2)	None	Good optical characterization scheme will be needed.
Aerogel radiator	K/pi separation between 3 to 10 GeV/c, e/pi separation ~2 GeV/c	· ·	Limited number of producers (one in Japan and one in Russian)	Extensive knowledge in the INFN group (led by Marco Contalbrigo) at CLAS12.



mRICH Key Components - Sensors & Readout (per module)

Components	Functions	Specs/ Requirements	Risk	Mitigation
Photosensor (MCP-PMT)	resolution; good timing	Magnetic field tolerant and radiation hardness (the same as for dRICH). 3mm x 3mm (or less) pixel size	Cost, R&D timeline, a single manufacturer for high quality volume production	R&D test needs to be performed in close collaboration Junqi Xie at Argonne, Mickey Chiu and Alexander Kiselev at BNL
Photosensor (SiPM - matrix)	Single photon detection and good spatial resolution (~ns)	Magnetic field tolerant and radiation hardness (the same as for dRICH). 3mm x 3mm (or less) pixel size		Proper cooling (tested in the 2nd mRICH beam test in 2018); further R&D is needed.
Readout	Amenica and change single	Low noise, time resolution		Closely working with JLab,
electronics	Amplify and shape single photon analog signal.	is ~0.5 ns; high density connectors; modular.	Chip development	INFN and Hawaii groups for readout development.

mRICH



Pros

Sweet momentum coverage for K/pi separation from 3 GeV to close to 10 GeV. It also provides the capability of e/pi separation around 2 GeV.

Modular design for array installation. Each module is independent with other modules and can be calibrated separately. Projective capability.

Performed two beam tests. The working principles have been validated in the first beam test in 2016 and the results have been published in NIM A. Further beam tests with tracking capabilities are expected and under planning.

Full GEANT4 simulation has been developed and verified using the beam test data.

An array of mRICH modules have been implemented in the sPHENIX for EIC simulations.

Cons

Photon sensors and readout electronics see direct hits of particles. Radiation hardness concerns.

Acquire aerogel tiles and maintain their long-term stability (optical)

Need high density photon sensors working in magnetic field.

Could create extra dead areas between modules. [Could be minimized by projective and creative integration schemes]

mRICH



Pros

Sweet momentum coverage for K/pi separation from 3 GeV to close to 10 GeV. It also provides the capability of e/pi separation around 2 GeV.

Modular design for array installation. Each module is independent with other modules and can be calibrated separately. Projective capability.

Performed two beam tests. The working principles have been validated in the first beam test in 2016 and the results have been published in NIM A. Further beam tests with tracking capabilities are expected and under planning.

Full GEANT4 simulation has been developed and verified using the beam test data.

An array of mRICH modules have been implemented in the sPHENIX for EIC simulations.

Cons

Photon sensors and readout electronics see direct hits of particles. Radiation hardness concerns.

Acquire aerogel tiles and maintain their long-term stability (optical)

Need high density photon sensors working in magnetic field.

Could create extra dead areas between modules. [Could be minimized by projective and creative integration schemes]

mRICH



Pros

Sweet momentum coverage for K/pi separation from 3 GeV to close to 10 GeV. It also provides the capability of e/pi separation around 2 GeV.

Modular design for array installation. Each module is independent with other modules and can be calibrated separately. Projective capability.

Performed two beam tests. The working principles have been validated in the first beam test in 2016 and the results have been published in NIM A. Further beam tests with tracking capabilities are expected and under planning.

Full GEANT4 simulation has been developed and verified using the beam test data.

An array of mRICH modules have been implemented in the sPHENIX for EIC simulations.

Cons

Photon sensors and readout electronics see direct hits of particles. Radiation hardness concerns.


Acquire aerogel tiles and maintain their long-term stability (optical)


Need high density photon sensors working in magnetic field.

Could create extra dead areas between modules. [Could be minimized by projective and creative integration schemes]

From a single module to large arrays












When particle is incident perpendicular at the center of mRICH, one can obtain the ring radius and the number of Cherenkov photons analytically. The details can be found in the appendix of the published mRICH NIM paper.

$$r = f \cdot \tan \theta$$
$$= f \cdot \sqrt{\frac{n^2 \beta^2 - 1}{1 - (n^2 - 1)\beta^2}}.$$

$$\begin{split} N_{\gamma} &= 2\pi\alpha L \left(1 - \frac{1}{\beta^2 n^2}\right) \\ &\times \int_{\lambda_1}^{\lambda_2} QE(\lambda) \cdot T_{aerogel}(\lambda) \cdot T_{lens}(\lambda) \cdot T_{glass\ window}(\lambda) \frac{d\lambda}{\lambda^2}. \end{split}$$







Modular focusing ring imaging Cherenkov detector for electron—ion collider experiments*


C.P. Wong g,*, M. Alfred i, L. Allison o, M. Awadi i, B. Azmoun c, F. Barbosa m, L. Barion j,r,


Examples from Analytical Calculations (perfect tracking!)












mRICH PID Parameterization




• We have implemented the pid performance (for particles entering mRICH perpendicularly as a starter) parameterization following Tom Hemmic's virtual PID base class: <mRICH.h>. This will be regised with more accurate results from GEANT4 simulations.

```
#ifndef __MRICH_H__
#define __MRICH_H__
// Created on 3/15/2020 by Xiaochun He at Georgia State University
// for the EIC Yellow Report
// Wrapper class for mRICH (fast PID for the EIC mRICH)
#include "PID.h"
#include "mRICHPidFast.cxx"
class mRICH: public PID
 mRICH(double trackResolution=0.5, double timePrecision=1.0, double pLow=3.0, double
 virtual ~mRICH() {}
 bool valid (TVector3 p) {return (p.Mag() > pLow && p.Mag() <
 double numSigma(TVector3 p, PID::type PID);
                (double numSigma, PID::type PID);
                 (double numSigma, PID::type PID) {r
                () {return myName;
  void description ();
  std::string myName;
  mRICHPidFast pid;
 mRICHPidInfo info;
 double fTrackResolution; // solution of the traker [mrad]
                          // time precision of the MCP-PMT [ns]
 double fTimePrecision;
 double pLow:
 double pHigh;
 double fSensorQMefficiency; // photon sensor quantum efficiency
 int mRICH_ID; //
```

mRICH Performance by set in GEANT4 simulation



- Projected e/pi separation of mRICH 2nd prototype detector (blue solid line)
- 2nd prototype detector can achieve 3sigma e/pi separation up to 2 GeV/c



- Projected K/pi separation of mRICH 2nd prototype detector (Green dots)
- 2nd prototype detector can achieve 3sigma K/pi separation up to 8 GeV/c



Backup slides

Resolution, Resolution, and Resolution



From
$$p = \gamma m v$$
 one gets, $m = p/(c\beta\gamma)$ and $\left(\frac{dm}{m}\right)^2 = \left(\frac{dp}{p}\right)^2 + \left(\gamma^2 \frac{d\beta}{\beta}\right)^2$.

In most cases, since γ is large, the mass resolution is determined mainly by the accuracy of the velocity measurement. The velocity measurement is given by

$$\frac{o_{\beta}}{\beta} = \tan(\theta_C) \sigma_{\theta_C}$$

With the average angular resolution for the angle photon σ_{θ_i} , the total resolution becomes

$$\sigma_{ heta_{C}}^{2} = \left(rac{\sigma_{ heta_{i}}}{\sqrt{N_{pe}}}
ight)^{2} + \sigma_{ heta_{Glob}}^{2}$$

The term $\sigma_{ heta_{Glob}}$ combines all contributions that are independent of the single photon measurement.

Resolution, Resolution, and Resolution



From
$$p = \gamma m v$$
 one gets, $m = p/(c\beta\gamma)$ and $\left(\frac{dm}{m}\right)^2 = \left(\frac{dp}{p}\right)^2 + \left(\gamma^2 \frac{d\beta}{\beta}\right)^2$.

In most cases, since γ is large, the mass resolution is determined mainly by the accuracy of the velocity measurement. The velocity measurement is given by

$$\frac{o_{\beta}}{\beta} = \tan(\theta_C)\sigma_{\theta_C}$$

With the average angular resolution for the angle photon σ_{θ_i} , the total resolution becomes

$$\sigma_{ heta_C}^2 = \left(rac{\sigma_{ heta_i}}{\sqrt{N_{pe}}}
ight)^2 + \sigma_{ heta_{Glob}}^2$$

 $\sigma_{\theta_c}^2 = \boxed{\frac{\sigma_{\theta_i}}{\sqrt{N_{pe}}}} + \sigma_{\theta_{Glob}}^2 \qquad \begin{array}{l} \text{The 2nd term is similar to Tom's equation, which} \\ \text{depends the external effects (background hits, tracking Parameters, etc).} \end{array}$

The term $\sigma_{ heta_{Glob}}$ combines all contributions that are independent of the single photon measurement.





$$\sigma_{\theta_i} = \sqrt{\sigma_{EP}^2 + \sigma_{Chro}^2 + \sigma_{Det}^2}$$


In the lens based mRICH design, σ_{EP} is minimized at the lens focal plane, σ_{Det} can be controlled with the lens focal length, and σ_{Chro} is reduced by selecting the lens transmittance in the near-UV region. As a consequence, the lens insertion improves the RICH performance at high momenta even in a compact device. Preliminary studies show that a 3σ kaon–pion separation is achievable at 10 GeV/c momentum with an aerogel of refractive index n = 1.03, a focal length of 6" and a detector pixel of 2 mm.



Tracking Requirement - Follow Tom's Example

- ▶ At any $|\vec{p}|$, there exists a "1-sigma" resolution in $\delta\theta_C$ which we'll call $\delta\theta_{C1}$
- Physics will demand that our full uncertainty follows: $\frac{\delta\theta_{C1}}{N}$.
- ▶ If we take the detector performance with a perfect tracker to be $\delta\theta_{C0}$.
- Then our external requirement on the tracker system is:

$$\left(\frac{\delta\theta_{C1}}{N}\right)^2 \le (\delta\theta_{C0})^2 + (\delta\theta_{tr})^2$$








120 GeV/c proton

















No precision tracking was available. Beam size is ~6mm in radius.



(b) (c)







18

(c)



No precision tracking was available. Beam size is ~6mm in radius.














No precision tracking was available. Beam size is ~6mm in radius.



















No precisic tracking wa available. size is ~6m radius.








Num. Of Photons On Ring


R: $40.82 \pm 1.08 \text{ (mm)}$ θ_c : $253.88 \pm 6.38 \text{ (mRad)}$ $\sigma_{\theta_c} = 6.38 \text{ (mRad)}$ 0 = 6.38 (mRad) 0 = 6.38 (mRad) 0 = 6.38 (mRad) 0 = 6.38 (mRad)




No precisic tracking was available. I size is ~6m radius.

120









GEANT4 Simulati


H (mm)



























Position scans with 120 GeV/c proton beam







Position scans with 120 GeV/c proton beam



Ring image from proton beam at an angle (110)






11 degree tilt downward

120 GeV/c proton beam






No precision tracking was available. Beam size is ~6mm in radius.



(b) (c)













No precision tracking was available. Beam size is ~6mm in radius.





(C)












No precision tracking was available. Beam size is ~6mm in radius.







GEANT4 Simulation


5/1/2020













Aerogel + fresnel lens. NO tech risks!

3 - 10 GeV/c k/pi, <2 GeV/c e/pi

YES

YES

YES, mainly to quantify the performance

Well advanced GEANT4 simulation for the standalone mRICH



Following points will be addressed to the best of my knowledge

- Technology used: spell out clearly any risk associated, if any Aerogel + fresnel lens. NO tech risks!
- Momentum range covered: p versus theta and Nsigma vs. p 3 10 GeV/c k/pi, <2 GeV/c e/pi
- Robustness of the design (e.g. sensitivity to magnetic field) and has a prototype been built?
- Are the electronics considerations clear (channel count, data size, rate, background) YES
- Time needed to complete the R&D and available workforce YES, mainly to quantify the performance
- Status of Simulation and Reconstruction Well advanced GEANT4 simulation for the standalone mRICH