
Testbeam analysis update

Abhisek Sen

Calibrations

Reconstructed energy for hadrons (EMCAL + HCALIN + HCALOUT)



Run(s)/sPHENIX Calo	Run(s)/PbGI, 1200V	Run(s)/PbGI, 1100V	Beam	Event	Comment
2722			-16	50K	EMCal normal bias C1 = 1.3 psia, C2 = 1.5 psia
2723, 2724		2725	-16	50K	EMCal lower bias @ gain of 1.15E5. C1 = 1.3 psia, C2 = 1.5 psia

- ➤ Needed a run dependent calibration scheme.
- Calibrations can depend
 - Running conditions
 - Energy (Sampling fractions changes)

Methodology

- * Created a root minimizer to give best possible reconstructed energy.
- * Tower-to-tower calibrations: HCAL: Cosmic, EMCAL: MIPs
- * Overall scale:

$$E_{reco} = p_1 E_{EMCAL} + p_2 E_{HCALIN} + p_3 E_{HCALOUT}$$
$$p = Min(\sum_{i=0}^{N_{events}} (E_{reco} - E_{truth}))$$


Used Minuit2:

ROOT::Math::Minimizer *min =

ROOT::Math::Factory::CreateMinimizer(minName, algoName);


- * Two steps:
 - EMCAL_MIP events: $E_{reco} = p_1 E_{HCALIN} + p_2 E_{HCALOUT}$
 - Total: $E_{reco} = p_3 E_{EMCAL} + p_4 (p_1 E_{HCALIN} + p_2 E_{HCALOUT})$

EMCAL MIP events




- Some improvement in the overall resolution.
- ❖ E_{EMCAL}<0.5GeV</p>
- Cherenkov cut < 10 to select a good hadron sample.</p>

EMCAL MIP events





- Same as before but now without the muon MIP peak.
- Overall the algorithm works better without the muon MIP peaks.

HCAL weight parameters




Resolution: EMCAL MIP







All three calorimeters



 $E_{reco} = p_3 E_{EMCAL} + p_4 (p_1 E_{HCALIN} + p_2 E_{HCALOUT})$

Resolution: ALL






Running conditions changed for last two energy runs.


Run(s)/sPHENIX Calo	Run(s)/PbGI, 1200V	Run(s)/PbGI, 1100V	Beam	Event	Comment
2692	2693		-2	>100K	EMCal normal bias. C1 = 12.0 psia, C2 = 12.5 psia
2695			-4	>200K	EMCal normal bias C1 = 12.0 psia, C2 = 12.5 psia
2696	2697		+4	>200K	EMCal normal bias C1 = 12.0 psia, C2 = 12.5 psia
2698 2700 2701	2702		-6	57205	EMCal normal bias C1 = 12.0 psia, C2 = 12.5 psia
2703,2704,2705,2706,2707,2708,2709,2710			-8	>200K	EMCal normal bias C1 = 5 psia, C2 = 6
2711,2712,2713,2714,2715	2716	2717, 2718	+8	>200K	EMCal normal bias C1 = 5 psia, C2 = 6
2719,2720,2721		2726	-12	>100K	EMCal normal bias C1 = 1.3 psia, C2 = 1.5 psia
2722			-16	50K	EMCal normal bias C1 = 1.3 psia, C2 = 1.5 psia
2723, 2724		2725	-16	50K	EMCal lower bias @ gain of 1.15E5. C1 = 1.3 psia, C2 = 1.5 psia
2727		2728, 2729	-24	>100K	EMCal lower bias @ gain of 1.15E5 C1 = 0.4 psia, C2 = 0.5 psia,
2730		2733, 2735, 2736, 2737	-28	>100K	EMCal lower bias @ gain of 1.15E5 C1 = 0.4 psia, C2 = 0.5 psia,


All three calorimeters



Same as before but without the muon peaks.

Resolution: ALL







Better resolution when relative calibrations are done without the muon MIP peak

Tilting


Normal position




HORIZONTAL POSITION




PRELIMINARY





+5 degree : EMCAL MIP



-5 degree : EMCAL MIP









Before calibration After calibration

+5 degree : ALL Calorimeters





-5 degree : ALL Calorimeter

