

# Testing MCPs in High $\vec{B}$ Facility

Kijun Park <sup>1</sup>,


C. Hyde<sup>1</sup>, G. Kalicy<sup>1</sup>, L. Allison<sup>1</sup>,  
Y. Ilieva<sup>2</sup>, T. Cao<sup>2</sup>, E. Brin<sup>2</sup>, C. Barbar<sup>2</sup>,  
P. Nadel-Turonski<sup>3</sup>, C. Zorn<sup>3</sup>

<sup>1</sup>Old Dominion University, <sup>2</sup>University of South Carolina, <sup>3</sup>Jefferson Lab

Oct.29,2015

# Update of MCPs testing at High $\vec{B}$ Lab


2013  
-2014



2015




# Magnet




- Superconducting solenoid
  - Diameter: 5-inch,  $L$ : 30-inch, warm bore
- $\vec{B}_{max} = 5.1\text{T}$  at 82.8A
  - Central field inhomogeneity:  $\leq 5 \times 10^{-5}$
  - Cylindrical volume ( $r=2.5\text{cm}, L=5\text{cm}$ )
- Non-magnetic stand for easy access to bore
- Mark down field map boundary on the ground
- Safety beacon's operated by magnetic field

# Magnet Cryogenics/Current Control




- 250 litter LHe Dewar allows a week measurement !  
→ Magenet (50l) is the driving mode
- LHe re-fill procedure  
→ Two persons work, they should communicate effectively  
→ Practical training is VERY important  
→ Follow the EXACT procedure and safety !  
→ Always monitor LHe level on magnet, in the Dewar
- Be familiar with accidental incident (We experienced couple of times)
- All detail procedure is on the logbook(both hard copy and Elog)
- Operator should be familiar with magnet current control

# Test Dark Box




- Non-magnetic, light-tight
- A cylindrical shape  
-  $D_{in} = 4.5\text{-inch}$ ,  $L=18\text{-inch}$ )
- Install a trail and turn-table allows to tranfort and rotate sensor
- LED light source (470nm) with pulser + optical fiber


# Test Dark Box



# Signal Processing




# Signal Processing






# MCP Measurement Strategy



Sweep magnetic field  
up to 1T with PMT210  
at certain angle set



- Test MCPs:  
→ PP0365g, PMT210/240, Katod, ...
- Minimize the Mechanical distortion  
→ MCP has a magnetizable component as part (coaval?) → Photonis(pp0365g):  
mechanical test up to 1T (signal disappear)  
→ Photek210: mechanical test up to 2T (signal disappear)
- Angle( $\theta$ ) scan with various increments
- Keep same  $\phi$  angle (reference position)
- Photek240 is too big to fit into a turn-table
- Pre-Amp( $\times 200$ ) used for both Photonis(pp0365g) and Photek210
- Feeding signal in fADC should be lower than 1V (jumper set: .5/1./2. V)



- login: Hall-B network  
→ ssh -X username@hallgw.jlab.org → ssh -X clasrun@highest2.jlab.org
- Start DAQ RunControl GUI : hbtest\_start/hbtest\_exit
- Connect → Configure (hbtest/hbtest\_er) → Download
- Prestart → Go

# DATA Storage

```
https://logbooks.jlab.org/entry/3292336
Photek PMT210 sensor measurement data table
Lognumber 3292336. Submitted by parkkj on Wed, 08/27/2014 - 09:58.
Last updated on Wed, 08/27/2014 - 18:06

Logbooks: HIGHB

*****
Sensor Configuration: B-FIELD TEST AT 30 DEG. ANGLE
Sensor: Photek PMT218
Angle(theta:along to z-axis)= 30 (deg)
HV set: -4300 V, Pre-Amp(x200)
*****

BIT: I B(A): RunNum: DataRate(Hz): EventNum(K)
0.0, 0.0, 241, 400, 347
0.1, 1.623, 240, 400, 151
0.2, 3.247, 249, 400, 162
0.3, 4.870, 251, 400, 151
0.4, 6.494, 252, 400, 151
0.5, 8.117, 253, 400, 218
0.6, 9.741, 254, 400, 175
0.7, 11.364, 255, 400, 153
0.8, 12.988, 256, 400, 153
0.9, 14.612, 257, 400, 151
1.0, 16.235, 258, 400, 150

*****
Above the 1.9T field signal is disappeared.

*****
Sensor Configuration: B-FIELD TEST AT ZERO ANGLE
Sensor: Photek PMT218
Angle(theta:along to z-axis)= 0 (deg)
HV set: -4300 V, Pre-Amp(x200)
*****

BIT: I B(A): RunNum: DataRate(Hz): EventNum(K)
0.0, 0.0, 259, 400, 205
0.1, 1.623, 260, 400, 169
0.2, 3.247, 261, 400, 232
0.3, 4.870, 262, 400, 175
0.4, 6.494, 266, 400, 161
0.5, 8.117, 264, 400, 130
0.6, 9.741, 268, 400, 235
0.7, 11.364, 271, 430, 226
0.8, 12.988, 273, 400, 198
0.9, 14.612, 274, 430, 223
1.0, 16.235, 275, 400, 218
1.1, 17.858, 276, 400, 231
1.2, 19.480, 280, 400, 187
```

```
https://logbooks.jlab.org/entry/3292340
Photonis PP0365G sensor measurement data table
Lognumber 3292340. Submitted by parkkj on Wed, 08/27/2014 - 10:22.

Logbooks: HIGHB

*****
Sensor Configuration: B-FIELD TEST AT ZERO ANGLE
Sensor: Photonis pp0365g
Angle(theta:along to z-axis)= 0 (deg)
HV set: -2350 V, Pre-Amp(x200)
*****

BIT: I B(A): RunNum: DataRate(Hz): EventNum(K)
0.0, 0.0, 323, 400, 183
0.1, 1.623, 325, 400, 201
0.2, 3.247, 326, 400, 226
0.3, 4.870, 327, 430, 155
0.4, 6.494, 328, 400, 155
0.5, 8.117, 329, 400, 162
0.6, 9.741, 330, 400, 154
0.7, 11.364, 332, 430, 151
0.8, 12.988, 333, 460, 184
0.9, 14.612, 336, 430, 205
1.0, 16.235, 340, 400, 307
1.1, 17.858, 343, 400, 161
1.2, 19.480, 342, 400, 153
1.3, 21.103, 343, 400, 166
1.4, 22.726, 344, 400, 156


*****
Above the 1.9T field signal is almost not visible and
Mechanical test done up to 2.6T

*****
Sensor Configuration: B-FIELD TEST AT 30 DEG. ANGLE
Sensor: Photonis pp0365g
Angle(theta:along to z-axis)= 30 (deg)
HV set: -2350 V, Pre-Amp(x200)
*****

BIT: I B(A): RunNum: DataRate(Hz): EventNum(K)
0.0, 0.0, 345, 430, 238
0.1, 1.623, 352, 430, 138
0.2, 3.247, 349, 430, 151
0.3, 4.870, 348, 430, 134
0.4, 6.494, 349, 430, 152
0.5, 8.117, 351, 430, 145
0.6, 9.741, 354, 430, 195
0.7, 11.364, 354, 430, 294
0.8, 12.988, 355, 430, 294
```


- DATA tables are available:
  - Photek: <https://logbooks.jlab.org/entry/3292336>
  - Photonis: <https://logbooks.jlab.org/entry/3292340>
- local computer : highest2.jlab.org
- ifarm (volatile disk) : /volatile/eic/commissioning/
- Tape Silo: /mss/home/parkkj/highBtest/

# Photek210 Drawing





- Single, Chevron, Gain  $\sim 10^6$  (2MCPs)
- FWHM 100ps

# Photek210 Voltage Divider




- Allow to provide voltage separately to Cathod-MCP, MCP-MCP, MCP-Anode
- Allow to combination of different voltage sets

# Selected Preliminary Result




# Selected Preliminary Result



- Photek210,  $3\mu\text{m}$  pore, Voltage:  $-4.76\text{kV}$ 
    - Normalized to the nominal point at 0T and 0(deg)
- Remark1)  $\theta = 0$  (deg),  
→ linearly decrease signal as increasing  $\vec{B}$
- Remark2)  $\theta > 0$  (deg),  
→ exponentially decrease signal  $\vec{B} > 0.5\text{T}$
- Remark3) No significant  $\phi$  angle dependent was observed


# Selected Preliminary Result



- Photek210, Modification voltage divider
- $\theta = 0$  (deg)
  - Change voltages: Cathod-MCP, MCP-MCP, MCP-Anode
  - Study of volatge correlation w.r.t gain
  - testing all possible combination voltages with various  $\theta$ ,  $\phi$ ,  $\vec{B}$  sets - MCP-MCP volatage change is the most sensitive
  - Overall preliminary uncertainty 5%
  - Recover gain at standard HV set by maximize HV(MCP-MCP)
  - Gain can be recoverable up to 5T




# Selected Preliminary Result



- Photek210, Modification voltage divider
- $\theta = 10$  (deg)
  - Control voltages: Cathod-MCP, MCP-MCP, MCP-Anode
  - Study of volatge correlation w.r.t gain
  - testing all possible combination voltages with various  $\theta$ ,  $\phi$ ,  $\vec{B}$  sets - MCP-MCP volatage change is the most sensitive
  - Overall preliminary uncertainty 5%
  - Recover gain at standard HV set by maximize HV(MCP-MCP)
  - Gain can be recoverable up to 3T

[\*\* Note: The plot is not normalized by  $\theta = 0(\text{deg})$  and 0T]

# Selected Preliminary Result



- Photek210, Modification voltage divider
- $\theta = 40$  (deg)
  - Control voltages: Cathod-MCP, MCP-MCP, MCP-Anode
  - Study of volatge correlation w.r.t gain
  - testing all possible combination voltages with various angle,  $\vec{B}$  sets - MCP-MCP volatage change is the most sensitive
  - Overall preliminary uncertainty 5%
  - Recover gain at standard HV set by maximize HV(MCP-MCP)
  - Gain can be recoverable up to 1T

[\*\* Note: The plot is not normalized by  $\theta = 0(\text{deg})$  and 0T]

# Summary

- The high  $\vec{B}$  facility has been established and fully functional
  - Superconducting magnet is working very smooth with various modes/fields
  - Own DAQ, Data process allows to have quick feedback of measurement
  - All data(raw/evio/root) have been stored JLab tape silo
- Data collections: Photek PMT210, 240, Photonis,...
  - Data contains: fADC, QDC, TDC
  - Converting raw data to ROOT/Ntuple/XML for analysis

- Photek PMT210: overall good response  $\theta = 0(\text{deg})$  up to  $\vec{B} = 5\text{T}$ 
  - Total collected charge is maximum at  $\vec{B} \sim 0.5\text{T}$
  - The collected charge shows exponential decrease as  $\vec{B}$  increase for  $\theta > 5(\text{deg})$
  - $\theta > 5(\text{deg})$ , good response up to  $\vec{B} = 2\text{T}$
  - The collected charge drops significantly between  $\theta = 0(\text{deg})$  and  $\theta = 5(\text{deg})$
  - $\theta > 5(\text{deg})$ , it shows a slow change
- Photek PMT240: overall good response  $\theta = 0(\text{deg})$  up to  $\vec{B} = 2\text{T}$

# Summary

- Increasing the potential difference across the channel plates is a powerful means to recover the loss in gain due to the effect of the magnetic field.
- Gain recovery is strongly correlated with the angle between the MCP and field axes:
  - The larger  $\theta$  angle, the more limited is the range of fields where the sensor can be operated at the same gate.
- Overall, other optimizations for gain recovery need to be implemented if the orientation of the sensor relative to the field varies significantly
- Data at 20 deg and -20 deg are under analysis.



# BACKUP

# Charge Calibration using Phillips7120



- Oscilloscope set(1GHz): sample integration : to calculate charge
- Set one of charge (Q) modes (300 pC) in Phillips7120
- ADJ. is set certain value (to make sure you can read -DC[V])
- Remove the horizontal offset by constant fit
- $Q = \frac{1}{R} \int V(t)dt$ , where  $R = 50\Omega$
- Calculation : use the "KaleidaGraph"
- Choose range before pulse (in order to remove offset)
- Fit it with constant take "mean" value
- Subtraction mean to pulse value  
Formula Entry:  $C_2 = C_1 - \text{mean}$
- Integration (macro) in time window (timing range is given by screen set)
- Pulse measurements done with 4 sets of voltage change
- $Q = 31.08 \text{ (pC/V)} \times V(V) + 1.597 \text{ (pC)}$

# Charge Calibration using Phillips7120



•  $Q \text{ [pC]} = (1.9 \pm 1.2) + (0.01883 \pm 0.00039) \times ADC_{ch}$



# E-LOGBOOK

https://logbooks.jlab.org/book/highb

JEFFERSON LAB ELECTRONIC LOGBOOK

Logged in as parkj (Logout)

Add content Logbooks Tags Useful Links Preferences Help/About

Search

**HIGHB**

**Ongoing Downtime Events**

Accelerator: Crew Chief On Call: Mike Aiken 630-7050 [188 hours 40 minutes]

**DAY Wednesday (27-Aug-2014)**

Lognumber	Date	Author	Title
3292340	10:22	parkj	Photonis PP0365G sensor measurement data table
3292336	09:58	parkj	Protek PMT210 sensor measurement data table

**DAY Thursday (21-Aug-2014)**

Lognumber	Date	Author	Title
3291811	13:46	parkj	PLANACON tube setup

**DAY Wednesday (20-Aug-2014)**

Lognumber	Date	Author	Title
3291885	10:36	parkj	LHe return to target group

**SWING Sunday (17-Aug-2014)**

Lognumber	Date	Author	Title
3291470	15:26	parkj	ProST magnet LHe level = 35cm

**OWL Saturday (16-Aug-2014)**

Date Picker: August 2014

Date Range  
Display Settings  
Tag Filters  
Book Filters  
Autorefresh

link :  
<https://logbooks.jlab.org/book/highb>