


TS3 Uncertainty


01-02-2016

Quantum Efficiency

The composition of a CCD is predominantly silicon, and therefore the sensitivity of the sensor to different wavelengths of light is dependent on silicon's inherent properties.





TS3



Uncertainty Measurements


- Repeatability
- Reproducibility
- Instrument Bias

Repeatability



Photodiode ratio measurements in cleanroom agree; FIL measurement is under investigation.

Reproducibility




Temp: -95°C


Gain: fe55 method

Instrument Bias: Lamp Drift

Time versus Current
FIL: 12 hour at 500nm



Time versus Current
FIL: 1min at 500nm



To Do

Repeatability:

- Measurements taken over a short time to capture the precision of the gauges.
- Measurements taken over days.
- Measurements taken over runs separated by months.

Reproducibility:

- TS3-2 and TS3-3 compared to each other.
- TS3-2 and TS3-3 compared to FIL stand.
- TS3-2 and TS3-3 compared to Vendor measurements (Absolute).

Instrument Bias:

Lamp and Housing, Off-axis parabolic mirror, Iris shutter, Filter wheel and filters, Integrating sphere, Dark space with baffles, Glass Dewar window.

Stray Light

Measurement method 1:

Dark cleanroom versus normal lighting.

Measurement method 2:

- 1) Measure the flux with and without a glass plate in the beam that blocks all radiation below 320 nm.
- 2) Compare the signal at 210 nm without the glass to the signal at 210 nm with the glass.
- 3) When the entire desired signal is blocked by insertion of a glass plate, what remains is scattering radiation.
- 4) This comparison is based on a taxing but real measurement scenario. The test is similar in principal to ASTM E387, Standard Test Method for Estimating Stray Radiant Power Ratio of Dispersive Spectrophotometers.

Measurement method 3:

Double Monochromators: The intensity of the light of other colors in the exit beam is referred to as the stray light level and is the most critical specification of a monochromator. The Double-Mono helps quantify this.