Beam Instrumentation for the JUNO Detector Calibration System

K. Heeger Yale University

Technical Considerations & Requirements


Beam transport pipe

- monitor position of beam to tune and optimize transport
- standard technology
- can be moved in and out of beam path
- cabling and readout outside beam pipe on deck above detector

Beam injection pipe


- custom-designed beam pipe without steering elements, possibly segmented beam pipe
- all detectors, readout, and cables need to be contained within beam pipe
- require detectors to be moveable out of beam
- desire detectors to be moveable along beam pipe

Beam Transport & Injection Pipe




Ref: McKeown

Injection Pipe



Ref: McKeown

Beam Monitors in Injection Pipe



Beam Transport Monitoring

Ionization Strip Chambers






Figure 1. Schematic view of the ionization strip c



ionization strip chambers



Beam Injection Pipe

Segmented Silicon Detectors

multi-pixel silicon p-ion-diode array for detection of beta electrons




Figure 7: FPD dartboard pixel pattern, surrounded by a guard ring and bias ring. Open circles show the points where electrical contacts are made.

Ref: 1404.2925

example: KATRIN focal plane detector manufactured by Canberra Belgium




Figure 8: Mounted detector wafer (back) with spring-loaded pogo pins and feedthrough flange. The inset is a close-up view of the detector compression scheme.

R&D Required

- identify monitor positions
 - required positions along transport pipe
 - desired optimum positions along injection pipe
- identify and adapt strip ionization chambers for transport pipe
- determine if injection pipe can be continuous or needs to be segmented
- determine required segmentation of Si detector in injection pipe (simulations and geometry)
- determine if detector needs to be movable along injection pipe
- develop customized readout package that can be integrated in beam pipe
- develop mechanical mechanism for "rotating" detector in and out of beam inside injection pipe
- test Si detector with electron gun