

# Coherent electron Cooling

## Proof of Principle

Instrumentation

Toby Miller

MAC Review  
8-10 December 2014



# Outline

## Introduction

Layout, Detailed overviews  
Electron Beam Parameters

## Electron Beam Instrumentation

### Sub-system descriptions

Current Monitors  
Position Monitors  
Profile Monitors  
Emittance Measurement  
Energy Spread Measurement  
Loss Monitors

## FEL IR Instrumentation

Intensity measurement  
Transverse Profile measurement  
Spectral Analysis

## RHIC Instrumentation

Wall Current Monitor  
BPM's, Schottky

## Summary

## Installation location RHIC Sector 2




# CeC PoP Facility Layout

Same Layout,  
4 views

Just fits in without  
moving any RHIC  
magnets.

Cable runs on the  
order of 200-300'

## RHIC Sector 2 Interaction Region



**BROOKHAVEN**  
NATIONAL LABORATORY


**Office of Science**  
U.S. DEPARTMENT OF ENERGY

December 8-10, 2014

# Electron Beam Transport Overview

## RHIC Sector 2 IR

RHIC 40GeV/n AU and 22MeV electron  
Merger location



**Goal of Diagnostics:**  
**Electron trajectory alignment for lasing**

**“Proof of principle”**  
**Minimal but sufficient diagnostics provided**



**BROOKHAVEN**  
NATIONAL LABORATORY

**Office of  
Science**  
U.S. DEPARTMENT OF ENERGY


December 8-10, 2014

# Layout for Phases 1 & 2

- **Phase 1: Test of the Gun & SRF cavity, measuring dark current with**
  - Profile Monitor
  - BLMs
- **Phase 2: Test with cathode & laser to produce beam, measuring with**
  - BPMs
  - Profile Monitors
  - BLMs
  - ICT (Current Monitor)
  - Faraday Cup
  - Emittance measured with multi-slit mask and Profile Monitor
- **Test independent of RHIC operations.**

## Instrumentation:

ICT = 1  
BPM = 2  
Profile Mon = 2  
Slit Mask = 1  
PMT BLM = 6  
Faraday Cup = 1




# Layout of Phase 2.1

- Addition of the 5-Cell and dog-leg injection line, repeat phase 1 measurements with higher energy beam.
- Test independent of RHIC operations.


## Instrumentation:

ICT = 2  
BPM = 4  
Profile Mon = 4  
Slit Mask = 1  
PMT BLM = 8  
PIN Diode BLM = 1  
Faraday Cup = 2





# Layout of Additions for Phase 3



# Electron Beam Parameters

## ELECTRON BEAM PARAMETERS

Energy	2 & 22 MeV
Charge per bunch	0.5 – 5 nC
Electrons per bunch	3 – 6x10 <sup>9</sup>
E-beam current (1nC)	78 $\mu$ A
Repetition rate	78 kHz
RMS Normalized Emittance	< 5 mm•mrad
RMS energy spread	<1x10 <sup>-3</sup>
RMS bunch length	10 psec
RMS transverse beam size	1 mm
E-beam power	1.7 kW

### Design and installation challenges due to SC gun and Linac:

- 10<sup>-11</sup> torr vacuum
- 200C bake
- Particulate processing, class 100 near SC cavities





# Instrumentation Parameters

Measurement		Required Resolution/ Accuracy	Instrument Resolution/ Accuracy	Method	Notes
2 to 22 MeV e-beam transport instrumentation parameters					
Electron Beam Position		30 $\mu$ m/ 100 $\mu$ m	10 $\mu$ m/ 50 $\mu$ m	BPMs (15mm buttons) Dual plane stations	Reuse Libera's from ERL, these have 703MHz BP filters, averaged data, not bunch-by-bunch
Cooling Section Beam Position	Electrons	30 $\mu$ m/ 100 $\mu$ m	10 $\mu$ m/ 50 $\mu$ m	BPMs with <b>SHARED</b> 15mm buttons Dual plane stations	BNL Zync Electronics VME boards ring <b>700MHz</b> front-end filter for averaging over many electron bunches
	Ions	100 $\mu$ m/ 100 $\mu$ m	10 $\mu$ m/ 50 $\mu$ m	BPMs with <b>SHARED</b> 15mm buttons Dual plane stations	BNL Zync Electronics VME boards ring <b>39MHz</b> front-end filter for averaging over many ion bunches
Transverse Profile		100 $\mu$ m (10% of 1 mm beam)	50 $\mu$ m	YAG screen & CCD <i>Measured at low power</i>	25mm YAG:Ce crystal with 2MP CCD GigE camera for 39px/mm resolution
Bunch Charge		5%/1%	1-10 pC noise floor	Integrating Current Transformer	Bergoz ICT, BCM-IHR. Expect bunch charge 0.5-5nC 0.1 – 7 $\mu$ s integration window @ 10kHz max rep rate
Beam Current		1 $\mu$ A	(TBD)	Dump Faraday Cup	Isolated dump as a Faraday Cup.
Beam Loss		10 $\mu$ A Loss limit	1 $\mu$ A	PMT detectors with JLAB VME electronics	Designed for a sensitivity of 0.1 – 6 $\mu$ A of beam loss and an integrated trip limit of 10 – 60 $\mu$ A•ms
Emittance		10%	10%	Plunging multi-slit mask & Profile Monitor	Two position mask, Horizontal & Vertical slits
Energy Spread		<1x10 <sup>-3</sup>	1% of max $\Delta p/p$	YAG in dispersive section	Based on 1mm beam size and 4X horizontal beam size growth under max $\Delta p/p$



# Outline

---

## Introduction

Layout, Detailed overviews  
Electron Beam Parameters

## Electron Beam Instrumentation

### Sub-system descriptions

Current Monitors  
Position Monitors  
Profile Monitors  
Emittance Measurement  
Energy Spread Measurement  
Loss Monitors

## FEL IR Instrumentation

Intensity measurement  
Transverse Profile measurement  
Spectral Analysis

## RHIC Instrumentation

Wall Current Monitor  
BPM's, Schottky

## Summary



December 8-10, 2014

# Integrating Current Transformer

For making periodic measurements of bunches or bunch trains from 0.1 $\mu$ s to 7 $\mu$ s long, located downstream of the gun, and at dump.

## **Bergoz ICT-CF6-60.4-070-05:1-H-UHV-THERMOE**

Integrating type, In-flange CT

Bergoz BCM-IHR electronics  
10kHz option  
Background Subtract  
Noise floor 1-10 pC  
Calibrated

### Mechanical details:


60.4mm ID  
40mm axial length  
Rad-Hard option  
Bakeable to 180C  
Separate bake-out zone  
Internal TC, type-E

### Requirements for bunch charge at 0.5-5nC:

5% accuracy  
1% resolution



## Signal processing timing diagram: Gate width <0.1 $\mu$ s up to >7 $\mu$ s



# Beam Position Monitors – Pick-ups & Electronics


Dual-plane button pick-ups

11 locations (2 at 2MeV, 9 at 22MeV)

15 mm diameter molybdenum button's

MPF Inc PN: A9111-2-CF

SMA Connector, SS housing & 4.5" CF flanges




BPM 15mm button




BPM Pick-up housing

## BNL Zync electronics VME board:





- VME Form Factor
  - Use RHIC Controls Infrastructure
- Configurable RF Section
  - Hadron or Electron Beam Options
- 4 x 400MSPS A/D Converters
  - 2 Planes of Measurement / Board
- Integrated Front End Computer
  - FEC & FPGA on Single Chip (Zynq)
- Ethernet Connectivity (x2)
  - Controls Network
  - High Speed Interface for Feedback

- Test results below at the ATF with buttons showed better than 100um accuracy and 10um precision



## Simulations show strong signals



**BROOKHAVEN**  
NATIONAL LABORATORY


**Office of Science**  
U.S. DEPARTMENT OF ENERGY

December 8-10, 2014


# Transverse Profile Monitor

## BNL Design

- YAG:Ce crystal, (100 $\mu$ m thick)
- Normal to beam; copper mirror behind crystal
- 25mm aperture
- Two position pneumatic plunger
- 2MP GigE camera (model “Manta” by AVT)
- 50mm fixed lens
- Tested to better than 50 $\mu$ m optical resolution
- Actuator separable from UHV linear motion feed through
- Laser focus assist
  - 405nm laser injected through dichroic beam splitter to excite YAG for focusing
  - 450nm Longpass edge filter protects CCD from backscattered 405nm light.
- Alignment laser for downstream multi-slit uses backside mirror on YAG holder.



YAG crystal & mirror holder (back-to-back copper mirrors)







# Emittance Slit Measurement

## - Transfer from ERL (need new mask)


- Low Power Operations Only
- Positioned 0.35 m upstream of profile monitor
- Tungsten Slit mask, optimized expected emittance:
  - 10 slits, 2mm separation, 200 $\mu$ m width in 1.5mm thick.
- New Dual axis design for Horizontal & Vertical measurement
- Replaces single axis mask originally installed.



### ANALYSIS:

An algorithm was developed for analyzing the image from a multi-slit mask for emittance measurement.

Future plans are to automate the image analysis for on-line processing and data logging.



Dual Station Actuator retrofitted for new dual axis mask.

Intensity Distribution at mask







Image on profile monitor after drift distance



# Energy Spread Measurements

- Max. Energy Spread:  $5 \times 10^{-4}$
- Beam Size (d): 1mm (dia.)
- Dispersion after *dipole* & *quad*
- $\sigma_{\delta-H}$  expected to be 4mm in dog-leg section
- $\Delta$  Beam Size  $_H$ : (4 – 1 mm) x 34.4px/mm  $\approx 100$ px
  - $\sigma_{\delta-H}$  spread over 100 pixels
  - ➔ Resolution: 1% of max  $\Delta p/p$

Image of 25.4mm YAG as projected onto CCD



- Camera resolution: 2MP
  - 1292 X 964 pixels
  - ( $px_h \times px_v$ )
- $Pitch_{YAG} = W_{CCD}/px_v$ 

$$= 964px/28mm$$

$$= 34.4px/mm$$

$$= 30\mu m/px$$






# Loss Monitors

## Beam Loss Monitors

- **PMT Detector** (similar to JLAB-CEBAF, FEL)
- (Hamamatsu R11558 PMT)
- Fast ( $\ll 1\mu\text{s}$ ), HV tunable for sensitivity adjustments
- 50dB DR (5nA – 100uA from tube)
- Limited coverage, small cross section

### BLM Electronics

- The new version of BLM electronics developed at JLAB
- Dual response:
  - Linear response for beam loss protection
  - Logarithmic response for instrumentation
- Two 8-ch modules installed support 14 PMT BLMs
- Designed for sensitivity to 0.1 – 6 $\mu\text{A}$  beam loss and 10 – 60  $\mu\text{A}\cdot\mu\text{s}$  integrated trip limit.
- Controls integration is underway.



### Features:

- Provide both machine protection and diagnostic functions.
- Instantaneous readback of beam loss.
- 16 bit digital output for integrating and logarithmic signals.
- Fast response,  $\ll 1\mu\text{s}$  response time for integrating, 10 nA<sup>2</sup> for log.
- Wide dynamic range ( $>50\text{ dB}$ ) for logarithmic signals.
- Built-in self test and onboard signal injection.
- FPGA controlled.
- Local data buffer for integrating and logarithmic signals.
- VME interface and fully integrated into EPICS.
- Pulse beam measurement and continuous monitoring.
- Low cost ( $\leq \$100$  per channel).

### Infrared Camera (FLIR A310)



Check for beam pipe heating, or losses other detectors

can't see.

Remote image display & storage, Ethernet communication.

New Beam Loss Monitor for 12GeV Upgrade  
J. Yan, K. Mhoney, ICALEPCS 2009a



**BROOKHAVEN**  
NATIONAL LABORATORY



December 8-10, 2014

# Outline

---

## Introduction

Layout, Detailed overviews  
Electron Beam Parameters

## Electron Beam Instrumentation

### Sub-system descriptions

Current Monitors  
Position Monitors  
Profile Monitors  
Emittance Measurement  
Energy Spread Measurement  
Loss Monitors

## FEL IR Instrumentation

Intensity measurement  
Transverse Profile measurement  
Spectral Analysis

## RHIC Instrumentation

Wall Current Monitor  
BPM's, Schottky

## Summary




December 8-10, 2014

# Infrared Diagnostics Layout for Wiggler Light

IR Diagnostics include:



- Power (over two ranges)
- Transverse profile
- Spectral Analysis



Simulations predict a 40mm radiation pattern at 15m from the FEL.

Wiggler Parameters	Value
Length (3 sections)	3 x 2.8 m
Period ( $\lambda_u$ ) / poles	4 cm / 3 x 60 + 2.5 (match)
Strength (K)	0.50
Wavelength	14 $\mu$ m
Optical Power	250 nW – 250 mW

Permanent magnet wigglers (three) built by the Budker Institute, 32mm gap, 0.134T peak field, 60 poles over principle length of 2.5m (+2.5 matching poles, over 30cm)



**BROOKHAVEN**  
NATIONAL LABORATORY


**Office of Science**  
U.S. DEPARTMENT OF ENERGY

December 8-10, 2014

# Wiggler Light Intensity Measurement

Change in the wiggler light intensity is a direct indication of FEL amplification, necessary for coherent cooling.

- **HBW Low Intensity IR Photodetector:**
  - Used to resolve variations from pulse to pulse in the 78kHz beam
  - Sensitivity as low as  $1\mu\text{W}$  over  $2 - 14\mu\text{m}$
  - BW up to 1MHz
  - TE cooled (HgCd)Te detector under GaAs lens
  - May employ lock-in amplifier to improve S/N ratio at 78kHz
- **Thermal Power sensor**
  - Used to measure average power in CW mode ( $T_{\text{resp}} = 1.8\text{s}$ )
  - $10\mu\text{W} - 3\text{W}$ , over 3 selectable ranges
  - Sensitivity over  $0.19 - 20\mu\text{m}$  through a 9.5mm aperture
  - 1.8 second response time




# Wiggler Light Transverse Profile

The profile of the FEL radiation beam is a good indication of the position of the electron beam and is expected to exhibit signs of when the electron and ion beams overlap. **Two methods of instrumentation:**



- **Pyrocam™ IV sensor array:**
  - Pyroelectric crystal integrating sensor
  - 100nm – 100μm sensitivity
  - 320 x 320 pixels on a 25 x 25 mm sensor
  - Responds to short pulses, uses a chopper in CW mode
- **Scanned Profile Measurement**
  - Low cost motorized stage for photodetector with pinhole mask
  - Fast response of the photodetector provides full scan in ~30sec.
  - Existing technique used at eLens for direct eBeam transverse profile scan

Pyrocam™ IV with  
pyroelectric sensor array



Sample image benefits from the 100 kilo-pixels over 1 in<sup>2</sup> sensor yields a spatial resolution of 12.8 pixels/mm

IR photo-conductive detector  
with 0.1mm Pinhole mask



Example image reconstruction & analysis from X-Y array of intensity data points (courtesy of RHIC Electron Lens Pinhole Scanning Profiler)



**BROOKHAVEN**  
NATIONAL LABORATORY

**Office of  
Science**  
U.S. DEPARTMENT OF ENERGY

December 8-10, 2014





# Wiggler Light Spectral Analysis

We expect to see variations in the spectral line width during the cooling process.

- Acton VM-504 Spectrometer:
  - Made by Princeton Instruments
  - 75-groove/mm grating, centered around  $14\mu\text{m}$  wavelength
- For spectral analysis over the transverse plane:
  - Princeton's imaging spectrograph, model SCT 320
  - Zero astigmatism yields high resolution

Acton VM-504 Spectrometer from Princeton Instruments.



# Outline

---

## Introduction

Layout, Detailed overviews  
Electron Beam Parameters

## Electron Beam Instrumentation

### Sub-system descriptions

Current Monitors  
Position Monitors  
Profile Monitors  
Emittance Measurement  
Energy Spread Measurement  
Loss Monitors

## FEL IR Instrumentation

Intensity measurement  
Transverse Profile measurement  
Spectral Analysis

## RHIC Instrumentation

Wall Current Monitor  
BPM' s, Schottky

## Summary



December 8-10, 2014



# RHIC Ion Diagnostics for CeC PoP

## 1. Align one ion beam through the Cooler

Use existing RHIC instrumentation routinely used for aligning colliding ion beams

- “DX BPMs” are located upstream and downstream of the CeC beam line.  
10’ s of microns resolution, average orbit
- Orbit stability feedback loops ensure beam position during the store.

## 2. Measure Cooling Progress


Use methods developed for successful Stochastic Cooling

Schottky Pick-ups

Wall Current Monitor (see next slide)

3KHz -6GHz BW

Anti-cooling: See when things go bad  
IPM’ s will see transverse effects  
Watch loss monitors carefully



Longitudinal Schottky spectra at 4 GHz before and after stochastic cooling. Blue is reference, red is narrowed spectrum after cooling. M. Brennan EPAC06 THPCH078





# RHIC Wall Current Monitor Simulation for CeC

## RHIC Wall Current Monitor Simulations (ideal conditions)


Au 40GeV/n, 1.5ns rms bunch length sigma,  $10^9$  ions/bunch

22MeV electron beam, 20ps rms, locked to center of ion bunch



Courtesy A. Fedotov

## Wall Current Monitor Simulated frequency response



0.05GHz

.5GHz

5GHz

Courtesy I. Pinayev



December 8-10, 2014

# Summary

---

CURRENTLY, all instrumentation for phase 1 (less loss monitors & Emittance Slit) is installed and online supporting the search for first beam.

OPEN ISSUES include:

- **Wiggler IR light diagnostics:**
  - profile & spectrum requirements are not yet defined
  - compounded by high cost solutions
- **BPM's:**
  - new electronics design with custom RF front-ends
  - ability to distinguish between ion & electrons on shared pick-up electrodes

LOOKING FORWARD, cables pulled for phases 2.1 & 3 are 90% complete and components for both phases are being assembled for installation along with loss monitors & Emittance slit for phase 2.1 in early 2015 and phase 3 in summer 2015.

We can apply experience gained at:

Energy Recovery LINAC for electron beam commissioning.

RHIC to help demonstrate coherent electron cooling for the first time.

