
### PHENIX Hadronic Calorimeter Tile Characterization and Testing

#### Jamie Nagle University of Colorado, Boulder




## CU Team Shawn Beckman Sebastian Seeds Theo Koblesky

#### Goals of Test Stand for Tile Characterization

- Optimize light collection via wavelength shifting fiber pattern and reflective coating
- Balancing of light response for:
  - (1) relative distance to nearest fiber
  - (2) change in sampling fraction as a function of depth
- QA testing of final production panels prior to installation in sPHENIX

#### Test Stand at University of Colorado


- 2-D Stepper Motor with Computer Control
- 2 SiPM readout with BNL designed preamp
- 1 PMT trigger
- Sr-90 source and cosmic capabilities
- Fast DAQ with full waveform recorded







#### **Test Configuration**





Test results indicate 15% more light within 0.5 cm of fiber location.

Note we had a relatively low light yield on the panels we have made ourselves.


As a basic test, last summer we machined a single scintillator with a groove for one straight fiber.

Yield Normalized to Mean Over Bins (4,257262 Counts), SiPM B






#### **Uniplast Manufactured Small Tile**



# Normalized Mean Voltage

#### Test Results

Pulse Averaged Test Scan



Similar 15% increase in light yield near fiber. Much higher light yield with Uniplast tile.

Location of PMT trigger 6

#### Full Size Tile Test Stand

We are building a second test stand capable of testing full 1 meter long tiles. All parts on hand – expected completion in 2 weeks.

Second test stand to use in Boulder on next set of Uniplast tiles. Eventually this test stand will move to Brookhaven for full QA tests.





#### <u>Simulation</u>

We have been in contact with CMS and other groups that have used embedded wavelength shifting fibers.

We are using suggested software to compare with experimental tests. We note that all groups say that one must test all configurations – not best option to rely on simulation packages.

http://gentitfx.fr/litrani/

There are now 2 versions of Litrani! This one, that we will call "the old Litrani", based on a simplified flat geometry, and a new one, called <u>SLitrani</u>, that is based on the <u>ROOT</u> geometry <u>TGeo</u>. This new version of Litrani is much more powerful since it imposes no restrictions to <u>TGeo</u>. Any setup which can be described by <u>TGeo</u> can be used in <u>SLitrani</u>. <u>SLitrani</u> is a new version, almost entirely rewritten from scratch. <u>CINT</u> macros are not compatible between the 2 versions.

The 2 versions will remain available, although no developments will be carried on the old Litrani. Have a look at the new <u>SLitrani</u>, and decides which one you intend to use!

#### <u>Summary</u>

- Test setup already evaluating existing panels
- Full 1 meter setup ready in 2 weeks for testing next full size panels
- Different fiber and reflective coating pattern comparisons can be quickly compared
- Test stand capable of full QA testing either in Boulder or at Brookhaven